

# Study Spotlight: Clinical outcomes of Mix & Match for ZEISS AT LISA<sup>®</sup> tri and ZEISS AT LARA<sup>®</sup>




Seeing beyond

## Methodology

Song et al. 2021



**Clinical outcomes of combined implantation of an EDoF IOL and a Trifocal IOL in a Korean Population.**



Monocular defocus curves of AT LARA 829MP and AT LISA tri 839MP and binocular mix & match curve

## Results

Journal of Ophthalmology (2021).

DOI: [10.1155/2021/9034258](https://doi.org/10.1155/2021/9034258)


- Binocular visual acuities at all distances were better in patients with combined IOL implantation than monocular visual acuities achieved with each IOL.
- Combined implantation of the AT LARA 829MP IOL and the AT LISA tri 839MP IOL showed similar visual acuity results at distance vision, and better results at an intermediate distance (-1.5D) and at near distance compared to monocular visual acuity results.

## Methodology

Mc Neely et al. 2021



**Visual outcomes and patient satisfaction 1 and 12 months after combined implantation of EDoF and Trifocal IOLs.**



Binocular defocus curve at 1 month.

## Results

Int Ophthalmol (2021).

DOI: [10.1007/s10792-021-01970-3](https://doi.org/10.1007/s10792-021-01970-3)


- Binocular defocus curve at 1 month was showing a peak in distance visual acuity at distance and a very good and stable visual acuity from 50 to 33 cm (-2.0 D to -3.0 D).
- Mean binocular UDVA at 1 month and 12 months was  $-0.08 \pm 0.07$  logMAR and  $0.09 \pm 0.06$  logMAR

## Methodology

Acar & Tabakci 2021



**Clinical outcome comparison: bilateral Trifocal vs. mix & match EDoF and Trifocal IOLs.**



Binocular defocus curves at 6 months.

## Results


Int Ophthalmol (2021).

DOI: [10.1007/s10792-021-01925-8](https://doi.org/10.1007/s10792-021-01925-8)


- No statistically significant differences in binocular UDVA and UNVA between both groups.
- Binocular defocus curve showed significantly better outcomes for AT LISA tri and AT LARA mix & match group from -2.0 D to -4.0 D.

## Methodology

Kim et al. 2022



**Comparison of visual outcomes after two types of mix & match implanted Trifocal EDoF and Trifocal IOLs.**



Binocular distance-corrected defocus curve. \* indicate significant difference

## Results

Graefes Arch Clin Exp Ophthalmol (2022).

DOI: [10.1007/s00417-022-05710-w](https://doi.org/10.1007/s00417-022-05710-w)

- In both groups, the binocular defocus curve showed a mean visual acuity of 0.2 logMAR or more between +1.00 D and -3.00 D of defocus.
- Between -1.00 D and -1.50 D of defocus (corresponds to intermediate vision) AT LARA-LISA group showed better visual acuity than the Triumph-HP group