

Transforming possibilities into realities.



ZEISS TIVATO 700


zeiss.com/tivato



Seeing beyond

Enhance precision and efficiency in the OR.

Today's microsurgical landscapes require advanced technology capable of performing complex procedures in simple, intuitive ways. Reliable performance for all of your patients is key.





For demanding surgical applications in ENT, Spine and P&R look no further than **TIVATO® 700** from **ZEISS**, the latest **Advanced Visualization System** from ZEISS. Understanding your daily workload, this fully integrated visualization system enhances usability and is built on technology ahead of its time.

During advanced surgical interventions and post-procedure workflows, **ZEISS TIVATO 700** is your intelligent partner that enhances all-around performance when it is needed most.

Achieve enhanced visualization during your workflow.


New treatment modalities, enhanced visualization, and brilliant broadcasting of procedures in the operating room (OR) -- the surgical landscape is changing at a rapid pace. Without the best tools and most intuitive technologies, your workflow can be disrupted inside and outside of the operating room.

Workflow-Enhancing Visualization enables you to have state-of-the-art modalities at your fingertips. ZEISS TIVATO 700 enables effortless use of advanced surgical visualization modalities and delivers superb image clarity for every procedure.


Visualization – enhanced.

Experience our fluorescence options.

Leveraging brilliant apochromatic optics from ZEISS, allowing you to venture into new fields including the use of **intra-operative fluorescence**¹ options. Assess patency of vessels joined by anastomosis with **ZEISS INFRARED 800**. Go even further with **ZEISS YELLOW 560** to visualize fluorescence-stained structures for additional clinical applications.



Vessel patency after anastomosis with ZEISS INFRARED 800²



Lumbar spinal tumor with ZEISS YELLOW 560³

Visualize with 4K.

Inside the OR, you can broadcast the procedure every step of the way in brilliant detail with the **integrated 4K camera technology**, while the large external monitor allows for easy viewing by your entire team.

¹ Please use the fluorescent agent as per the approval status for the application in your country.

² Image courtesy of Prof. Dr. Milomir Ninkovic, University Clinic, Bogenhausen, Germany; recorded with ZEISS OPMI PENTERO 900

³ Image courtesy of Prof. Dr. Talat Kiris from LIV Hospital Cancer Center, Istanbul, Turkey



Always in position to perform.

The OR is a demanding place full of activity – with multiple people and devices all inside a small room. Yet, different types of procedures require different set-ups and workflows. Every day in the OR is different and any disruption, such as elaborate rebalancing, can hinder desired clinical outcomes.


Ultimate Reach & Flexibility gives you and your team the freedom to handle every OR situation with confidence.

ZEISS TIVATO 700 offers **excellent overhead clearance** to work under the arm of the visualization system in every set-up. Position the base as your space allows and still operate comfortably with **Ultimate Reach**. When working with long instruments, rely on the **large working distance** of 200 – 625 mm to ensure continuous focus. Prepare the device with ease for daily use with convenient features like **AutoBalance** and **AutoDrape®**.

Flexibility – reached.

Stability reimagined.

Avoid unwanted vibrations during procedures after the repositioning of the visualization system and during small movements. Experience the innovative **Active Vibration Damping**, strengthened by robotic movements, to streamline your surgical workflow.







Get fully connected in your OR.

Frequently rotating OR teams require easy-to-use medical equipment with intuitive user interfaces. Integrated medical devices lead to a faster turnaround in the OR enabling more patients to be treated and improved procedural efficiency.

Go **All-Digital** with an easy-to-use and future-ready platform for a seamless entry into the world of fully connected visualization systems.

ZEISS TIVATO 700 offers you an **intuitive graphical user interface** that makes interaction with the system simple for everyone in your operating room. It comes with a **full connectivity package**, including the **ZEISS Connect App** to digitally manage surgical data, and dedicated functionalities to integrate into an existing hospital IT infrastructure. The **ZEISS Observe App** enables streaming in realtime on a mobile device and ZEISS VR ONE Plus glasses. Remote connectivity, enabled by **ZEISS Smart Services** offers you faster support to improve availability of your device.

Digital ORs – connected.



ZEISS Connect



ZEISS Observe

Succeeding. The smart way.


Smart Services from ZEISS



Be successful. With remote solutions from ZEISS.

ZEISS TIVATO 700 is powered by **ZEISS Smart Services**, the new remote service platform from ZEISS. With this platform ZEISS strives to establish a new dimension in services. Its architecture provides a high level of IT infrastructure security and offers all stakeholders more confidence in conducting daily business.

Is maintaining system availability of critical importance to you? If this is the case, you should opt for fast, qualified assistance with an OPTIME service agreement including ZEISS Smart Services.



Technical Data

TIVATO® 700 from ZEISS

		Packages				
		Essential	Standard	Spine Comfort	ENT Comfort	Premium
System	TIVATO 700 Essential system	●	-	-	-	-
	TIVATO 700 basic system	-	●	●	●	●
Main Tubes	Straight binocular tube	○	○	□*	□*	□*
	180° tiltable tube	○	○	-	○	-
	Foldable Tube f170/f260	-	-	●	○	●
Eyepieces	10x wide-field eyepieces	○	○	○	○	○
	12,5x wide-field eyepieces	○	○	○	○	○
Focus	SpeedFokus	-	□	●	□	●
Working Distance	200–500 mm	●	-	-	-	-
	200–625 mm	-	●	●	●	●
Illumination	Xenon – Fully integrated Xenon light source with automatic lamp exchange	-	○	○	○	●
	TriLED – Fully integrated LED light source	●	○	○	○	-
Face-to-face observation	Face-to-face with StereoBridge	□	□	●	□	□
Side-Observation	Left/right with stereo co-observer	○	○	○	●	○
	Side Port Preparation left/right	○	○	○	-	○
	No Side Port – side observation excluded	○	○	○	-	-
Balance	GuidedBalance	●	-	-	-	-
	AutoBalance	-	●	●	●	●
XY-Movement	Angular movement with joystick or foot pedal	-	□	□	●	●
Video System	Standard HD Video Package – Fully integrated 1-chip HD camera, 1080p	●	○	-	-	-
	Premium HD Video Package – Fully integrated 3-chip HD camera, 1080p	□ ¹	○	○	○	-
	4K Video Package – Fully integrated 3-chip 4K camera, 2160p	-	-	○	○	●
Video Options	Video recording	□	□	□	□	●
	HD Video-In	-	□	□	●	●
Fluorescence**	YELLOW 560	-	-	□	○	●
	YELLOW 560 LED	□	□	□	□	-
	INFRARED 800	-	-	□	○	●
Visualization	Additional Depth of Field	□	□	□	●	●
Networking & Storage	Integrated 1TB harddrive	□	□	□	□	●
	Wireless Network Package – Enables WLAN and WiFi Hot Spot functionality	●	●	●	●	●
	LAN Interface ²	-	-	□	□	□
	Shared Network Package – Storage on shared network directories	□	□	□	□	●
	DICOM Package – Patient data exchange with PACS system via DICOM	□	□	□	□	●
	ZEISS Smart Services Package – Connectivity for ZEISS Smart Services	●	●	●	●	●
Navigation	Navigation Package	-	-	□	□	□

* available as an additional tube

** separate module necessary for activation

¹ only available with YELLOW 560 LED

² only available with Navigation Package

● Basic ○ Configurable □ Option

Cochlea implantation surgery: Image courtesy of Prof. Dr. Joachim Hornung, University Clinic Erlangen, Germany (Page 2 & 8)

Spinal hemangioblastoma: Image courtesy of Prof. Dr. med. Karl-Michael Schebesch, University Clinic Regensburg, Germany (Page 5)

Decompression surgery: Image courtesy of Dr. Andreas Korge, Schoen Clinics Munich, Germany (Page 6)



Carl Zeiss Meditec AG

Goeschwitzer Strasse 51–52

07745 Jena

Germany

www.zeiss.com/med/contacts

www.zeiss.com/tivato

en-INT_30_010_01271 Printed in Germany. CZ-X/2023

The contents of the brochure may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information. Subject to changes in design and scope of delivery and due to ongoing technical development, TIVATO, AutoDrape, AutoBalance, Varioskop and VisionGuard are either trademarks or registered trademarks of Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and/or other countries.

© Carl Zeiss Meditec AG, 2023. All rights reserved.