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
Time	User	Category	Description
15/03/2018 16:11:31	Operator	Execute	In task 1 'Image Processing' the following parameter has changed: 'Contrast' has changed from 1.05 to 1.11.
15/03/2018 16:11:31	Operator	SettingsChange	In task 1 'Image Segmentation' the following parameter has changed: 'Image Segmentation' the following parameter has changed: 'Image Segmentation Intersect Segment ParameterThreshold Value' has changed from 120 to 126.
15/03/2018 16:11:47	Operator	Execute	In task 1 'Image Segmentation' the following parameter has changed: 'Image Segmentation Intersect Segment ParameterThreshold Value' has changed from 126 to 126.
15/03/2018 16:12:04	Operator	SettingsChange	In tool 'AnalysisTaskGroupBase' the following parameter has changed: 'Intersect Segment ParameterThreshold Value' has changed from 120 to 126.
15/03/2018 16:12:19	Operator	Execute	Beginning loop iteration 3.
15/03/2018 16:12:19	Operator	Execute	In task 1 'Load Image' the following parameter has changed: 'File Name' has changed from 'C:\Fentic\data\Fentic-Feafile--02.tif' to 'C:\Fentic\data\Fentic-Feafile--03.tif'.
15/03/2018 16:12:38	Operator	Execute	In task 1 'Image Processing' the following parameter has changed: 'Contrast' has changed from 1.11 to 1.11.
15/03/2018 16:12:57	Operator	SettingsChange	In task 1 'Image Segmentation' the following parameter has changed: 'Image Segmentation Intersect Segment ParameterThreshold Value' has changed from 120 to 126.
15/03/2018 16:13:09	Operator	SettingsChange	In task 1 'Image Segmentation' the following parameter has changed: 'Image Segmentation Intersect Segment ParameterThreshold Value' has changed from 126 to 126.
15/03/2018 16:13:21	Operator	Execute	Beginning loop iteration 4.
15/03/2018 16:13:28	Operator	Execute	In task 1 'Load Image' the following parameter has changed: 'File Name' has changed from 'C:\Fentic\data\Fentic-Feafile--03.tif' to 'C:\Fentic\data\Fentic-Feafile--01.tif'.
15/03/2018 16:13:33	Operator	SettingsChange	In task 1 'Image Segmentation' the following parameter has changed: 'Image Segmentation Intersect Segment ParameterThreshold Value' has changed from 120 to 126.
15/03/2018 16:13:50	Operator	Execute	In task 1 'Image Segmentation' the following parameter has changed: 'Image Segmentation Intersect Segment ParameterThreshold Value' has changed from 126 to 126.
15/03/2018 16:14:01	Operator	SettingsChange	In tool 'AnalysisTaskGroupBase' the following parameter has changed: 'Intersect Segment ParameterThreshold Value' has changed from 120 to 126.
15/03/2018 16:14:22	Operator	Execute	4 loop iterations were executed.
15/03/2018 16:16:54	Operator	SetStatus	The job result 'C:\Users\jrbarnet\AppData\Local\Carl Zeiss\ZENCore\Results\20180315_160239_GxPGrainIntercept_3.rpt' was saved and due to GxP also signed using certificate 'L:', 'C:', 'CN=Operator' (3104682802051F4E3778A2802838E6F6).

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JP_42_013_203 | CZ 11-2019 | 構成、仕様、外観等は予告なく変更する場合があります。 | © Carl Zeiss Microscopy GmbH
Carl Zeiss Microscopy GmbH