

Deepen Your Knowledge.



ZEISS Crossbeam Family

Your FIB-SEM for High Throughput 3D Analysis and Sample Preparation

zeiss.com/crossbeam



Seeing beyond

Your FIB-SEM for High Throughput 3D Analysis and Sample Preparation

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › Technology and Details
- › Service

ZEISS Crossbeam combines the powerful imaging and analytical performance of a field emission scanning electron microscope (FE-SEM) column with the superior processing ability of a next-generation focused ion beam (FIB).

Crossbeam gives your 3D work that dynamic edge, whether you are milling, imaging or performing 3D analytics. Extract true sample information from your SEM images using Gemini electron optics. The Ion-sculptor FIB column introduces an altogether new way of FIB-processing. By minimizing sample damage you'll maximize sample quality—and perform experiments faster at the same time.

Customize your instrument to achieve both high quality and high throughput TEM sample preparation. Exploit the variable pressure capabilities of Crossbeam 350. Use Crossbeam 550 to prepare your most demanding samples, choosing the chamber size that best suits your samples. Or choose Crossbeam 550 Samplefab for advanced, fully-automated, unsupervised multi-site TEM lamella preparation.

You may be working on your own or in a multi-user facility, as an academic or in an industrial lab. If you've set your sights on high impact results, Crossbeam's modular platform concept lets you upgrade your system as your needs grow.



Simpler. More Intelligent. More Integrated.

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Maximize Sample Insights in Both 2D and 3D

Count on excellent images from any sample thanks to the Gemini electron optics of your ZEISS Crossbeam. You will achieve high resolution and contrast while reaping the benefits of high signal-to-noise ratios, right down to very low accelerating voltages. Prepare high quality samples, like TEM lamellae, using the FIB's low voltage performance and characterize your samples comprehensively in 3D. Use a wide choice of detectors, including the unique Inlens EsB (energy selective backscatter) detector for pure material contrast. Investigate non-conductive specimens undisturbed by charging artifacts, offset either with local charge compensation while keeping high vacuum in the chamber or with variable pressure available in Crossbeam 350.

Increase Your Sample Throughput

Combine Gemini optics with a new way of FIB machining: the superior low voltage performance of the Ion-sculptor FIB column delivers fast and precise results while keeping amorphization damage on your sample to a minimum. Use these advantages especially for the preparation of TEM lamellae – even challenging samples. Benefit further from the FIB's high current capability that saves time and achieves excellent FIB profiles with up to 100 nA current—without compromising the ultimate FIB resolution. Save even more time with automatically prepared batches of cross-sections or with any user-defined pattern. And the benefits just keep on coming throughout your long-term experiments as optimized routines enhance FIB source lifetime and stability.

Experience Best 3D Resolution


Enjoy precise and reliable results in FIB-SEM tomography with best 3D resolution and leading isotropic voxel size. The Inlens EsB detector lets you probe and image less than 3 nm in depth. Expand the capacity of your Crossbeam with Atlas 5, our market-leading package for fast, precise tomography. You will save time by collecting your serial section images while milling. You also have the advantage of trackable voxel sizes and automated routines for active control of image quality. Meanwhile, Atlas 5's new integrated Analytics module enables 3D EDS and 3D EBSD analysis during tomography runs.



Count on excellent images thanks to Gemini electron optics.



Benefit from the superior low voltage performance of the Ion-sculptor FIB-column, especially for TEM lamella preparation.



▶ [Click here to view this video](#)

3D tomogram of a solid oxide fuel cell anode made of the heat resistant composite material Nickel Samaria-doped ceria.

Your Insight into the Technology Behind It


- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Profit from Gemini Optics


Crossbeam's FE-SEM column is based on Gemini electron optics. You will appreciate the long-term stability of your SEM alignment and the effortless way it adjusts all system parameters such as probe current and acceleration voltage. Unlike other FE-SEMs, Gemini optics don't expose your specimen to a magnetic field. This allows you to achieve distortion-free, high resolution imaging over large fields of view as well as to tilt the specimen without influencing the electron optical performance. Even magnetic samples can be imaged easily.

Choose between Two Columns:


- The Gemini VP column of Crossbeam 350 gives you maximum sample flexibility and multi-purpose environments. With the optional Variable Pressure (VP) you can perform *in situ* experiments under excellent analytical conditions, even with outgassing or charging samples.
- The Gemini II column of Crossbeam 550 has a double condenser system that enables high resolution, even at low voltage and high current. It's ideal for high resolution imaging at high beam current and for fast analytics
- Simultaneous Inlens SE and Inlens EsB imaging provides unique topographical and material contrast. That means you will gain more information in less time.



Magnetic field leakage of the Gemini lens compared to a traditional single-pole lens design. The minimum magnetic field on the sample allows highest ion and electron beam performance on a tilted sample as well as high resolution imaging of magnetic materials.



ZEISS Crossbeam 350: Gemini column with single condenser, two Inlens detectors and VP capability.




ZEISS Crossbeam 550: Gemini II column with double condenser and two Inlens detectors.

Your Insight into the Technology Behind It


- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Discover a New Way of FIB-Machining – From Massive Ablation to Nanometer Precision

Maximize sample quality by using the low voltage capabilities of the Ion-sculptor FIB column. Minimizing amorphization of delicate specimens will give you the best results after thinning or polishing—with the added advantage of fast probe current exchanges to accelerate your FIB applications. Or opt for high current performance and double the speed of your 3D FIB-SEM applications by working with the high gallium ion beam currents. You'll get precise and reproducible results with maximum stability during the acquisition time. The column design gives you access to five orders of magnitude in beam current, from 1 pA up to 100 nA. The larger beam currents of up to 100 nA allow fast and precise material removal and milling processes. Meanwhile, at low currents you will achieve exceptionally high FIB resolution of less than 3 nm. The gallium focused ion beam source – the so-called LMIS (liquid metal ion source) Ga source – is designed for a typical lifetime larger than 3000 μAh when operating at a target emission current of 2 μA . For long-term experiments, you have the bonus of the Crossbeam family's automatic FIB emission recovery.



Regulation characteristic of Ga source emission current. The emission is stable for more than 72 hours.



ZEISS Crossbeam 550: FIB- and FE-SEM column arranged at an inclination angle of 54°.


Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service



Customize Your Crossbeam with the Remote Application Programming Interface

Innovative experiments will often require new functionality beyond what is provided by the operating software of your electron microscope. That's why the open programming interface of Crossbeam is designed to allow access to almost every microscope parameter. The remote API lets you take complete control of electron and ion optics, stage, vacuum system, detectors, scanning and image acquisition from custom programs – whether running on the system PC or on a remote workstation.

ZEISS provides both documentation and code examples in various programming languages – plus technical support to make sure you get the results you want. Quickly.



Workflow for lathe milling, implemented in the custom application SmartLathe and using the API interface.



Pillar for compression testing after being machined using lathe milling: SEM top view (left), SEM side view (right).


Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

How a Laser FIB Workflow Accesses Deeply Buried Structures


To access deeply buried regions, you need to localize ROIs in 3D, ablate material via a targeted preparation and do 3D imaging and analytics. Add a femtosecond laser to your ZEISS Crossbeam and benefit from ultra-fast sample preparation.

- Gain rapid access to deeply buried structures
- Prepare extremely large cross-sections up to millimeters in width and depth
- Benefit from minimal damage and heat affected zones due to ultrashort laser pulses
- Perform laser work in a dedicated chamber with debris handling to avoid contamination of the main instrument
- Find your hidden ROIs by correlation with previously acquired X-ray microscopy datasets
- Automate laser processing, polishing, cleaning and transfer of the sample to the FIB chamber




[▶ Click here to view this video](#)


The LaserFIB enables you to optimize and automate laser processing.



Deep laser cut in electronics sample to gain access to buried ROI in 860 µm depth (left). The targeted structures were already visible after laser preparation. Minimal time was needed to FIB polish only the ROI area (right) and reveal finest details of the microbump.



200 µm wide cross section in a ceramic sample with 200 µm clearance on each side cut by fs-laser in less than 30 s. This pattern can be used to investigate the micro-structure of the sample material in cross section or as pre-preparation for a subsequent FIB-SEM tomography run.



Surface detail of a cross section (similar to that shown on the left) in a metal alloy sample depicting the quality of the cut after polishing only with the laser. No FIB polish was applied. The grain structure as well as inclusions are directly visible on a large area.

Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

It's Easy to Create TEM Lamellae

For conventional TEM sample preparation, simply use the module ASP (Automatic Sample Preparation): it includes all necessary steps for a guided, semi-automated workflow and delivers TEM lamellae ready for lift-out.

How it Works

- Click on the TEM sample icon.
- Draw a line to define the location of the lamella on your sample.
- Trigger execution.

Prepare Batches of TEM Lamellae

- Execute a batch of TEM samples at predefined sites without supervision.

How it Works

- Define location of single TEM lamella and transfer to process list.
- Repeat step 1 as often as required or perform copy & paste in sample mode.
- Execute process list.



Automated chunk milling of TEM lamella in bulk sample according to a user-defined recipe is available on every Crossbeam. The TEM grid sample holders and loading accessories make it easy to handle your valuable TEM lamellae.





Illustration of the simple three-step workflow for TEM sample preparation (depo stands for deposition, DC for drift correction)



Array of TEM lamellae prepared automatically.


Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Select Your Micromanipulator for TEM Lamella Lift-out

When starting a TEM lamella preparation workflow find your region of interest quickly with the help of the super-eucentric 6-axis stage and always stay at eucentricity when tilting the sample no matter which working distance. Prepare your sample and finally, utilize a micromanipulator for the next steps in the workflow.


It's quick and easy to lift out a prepared TEM lamella from the bulk. Select a micromanipulator that is targeted to your needs in flexibility, freedom of operation and ease-of-use in control. Attach your lamella to a grid for final thinning and low kV polishing.



Find your region of interest quickly with the help of the super-eucentric 6-axis stage and always stay at eucentricity when tilting the sample. Prepare your sample. Finally, lift out your TEM lamella.



The micromanipulator of your choice will be configured to enable optimized workflows.




Attach the needle of the micromanipulator to the lamella, lift it out of the bulk and attach it to the TEM-grid for further investigation in transmission mode. (from left to right)

Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service


Keep Control During TEM Lamella Thinning

The final polishing step is crucial, as it defines the quality of your TEM lamella. When you aim to reach a desired thickness, the SEM allows live monitoring of the thinning. During imaging, the 'split mode' gives you the benefit of having the signals of several detectors available in parallel. Use the SE signal to judge lamella thickness and obtain reproducible end thickness. At the same time, the InLens SE signal helps you control surface quality.



Keep full control during lamella thinning.

SmartEPD is an optional software module that allows you to determine the thickness of your TEM lamella quantitatively and thus stop the thinning process at your pre-defined endpoint, this time using the InLens EsB detector.



Use SmartEPD to determine the lamella thickness and endpoint of polishing quantitatively.

Even more benefits come your way with the X²-Holder, which is built to enable the preparation of ultra-thin lamellae. This is a big help when dealing with challenging samples that show intrinsic stress, for example, heterogeneous materials or polymers that would otherwise bend.



[Click here to view this video](#)

Use the patented X² preparation method to prepare ultra-thin, stable TEM lamellae and obtain a homogenous thickness of less than 10 nm without causing sample damage.

Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Fully automated TEM preparation with Crossbeam 550 Samplefab

Achieve fully-automated and unattended TEM sample preparation in your semiconductor lab. Obtain best sample quality and take advantage of high automation reliability and lamella success rates, especially for multi-site preparations. Crossbeam 550 Samplefab is a robust, pre-configured high-end FIB-SEM for the automation of single and multi-site TEM lamella preparation.


For more information refer to
www.zeiss.com/samplefab

What You Can Expect:


- Exceptional lamella quality using the excellent low-kV performance of the Ion-sculptor FIB column.
- Highly accurate endpoint and thickness control on even the smallest semiconductor devices enabled by uncompromised SEM live imaging while milling.
- Recipes that can be optimized to each individual sample and can be saved for regular use.
- An intuitive user interface designed for fast learning and optimum efficiency without sacrificing flexibility – providing truly hands-free automation for long multi-sample runs.

Automated Workflow


The three main steps of TEM lamella preparation – chunking, *in situ* lift-out, and thinning – are enabled by the automated workflow. Run each of these steps independently or linked and create user-defined workflows for single or multiple TEM lamellae from only one or even multiple samples on your sample holder. Linking the chunking and thinning main steps provides TEM lamella ready for *ex situ* lift-out.



Chunking: Fiducial marks define each lamella location. Each lamella site can be assigned its own optimized recipe depending on the needs of the sample and location.



Lift-out: Crossbeam 550 Samplefab provides automation for *in situ* and *ex situ* lift-out workflows. Automated *in situ* lift-out, when combined with automated chunking, offers a robust hands-free “bulk to TEM grid” workflow with no operator intervention.



Thinning: Auto-thinning down to 100 nm lamella thickness is possible, both in bulk (for *ex situ* lift-out workflows) and on the TEM grid following lift-out. Endpointing relative to operator-defined fiducials keeps the region of interest in target.

Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

ZEISS Atlas 5 – Master Your Multi-scale Challenge

Use Atlas 5's sample-centric correlative environment to create comprehensive multi-scale, multi-modal images. This powerful yet intuitive hardware and software package extends the capacity of your Crossbeam. With its efficient navigation you can correlate images from any source.

For example, use X-ray volume data from your ZEISS X-ray microscope to target buried features of interest and analyze them in your Crossbeam. Take full advantage of Atlas 5's high throughput and automated large area imaging. Unique workflows will help you gain a deeper understanding of your sample. The modular structure lets you tailor Atlas 5 to your everyday needs in materials or life sciences research.




Recommended Modules for Your Crossbeam

- NPVE Advanced (Advanced Nanopatterning & Visualization Engine): Perform nanopatterning with full control over patterning geometry and parameters.
- 3D Tomography: Turn your Crossbeam into a precise 3D FIB-SEM tomography acquisition engine with automated sample preparation. Automatically acquire 3D image data with up to several thousand images and a voxel resolution below 10 nm isotropic voxel size in 3D. Unique sample-tracking technology gives you the benefit of consistent slice thickness over long acquisitions. Meanwhile, robust autofocus and auto-stigmation algorithms keep all of your images sharp.
- Analytics: Add 3D EDS / 3D EBSD analytics to high resolution FIB-SEM tomography acquisition. Specify imaging and mapping conditions independently. Use the advanced acquisition engine to automatically switch between analysis and imaging conditions during acquisition. Flexible visualization allows you to simultaneously view SEM images and process elemental maps.



▶ [Click here to view this video](#)



▶ [Click here to view this video](#)

Tomography data of a lead free solder containing Cu and Ag particles in an Sn matrix. SEM images (top) and EDS maps (bottom) were acquired at the same sample site at 1.8 kV and 6 kV, respectively, with a ZEISS FIB-SEM and Atlas 5 Analytics. Courtesy of: M. Cantoni, EPFL Lausanne, Switzerland.

Expand Your Possibilities


- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Fastmill – Speed Up Your Material Removal

Milling speed depends on a multitude of factors: target material, lattice orientation, ion current, milling geometry and so on. For a given material, your scanning strategy has the biggest impact on the material removal rate. During milling, sample topography changes are based on the precise milling strategy. This change, in turn, affects the milling rate.

Two milling styles are commonly distinguished: line and frame milling. In the first, the ion dose is delivered line by line in a single pass. In the latter, the entire frame is milled multiple times until the total dose is delivered. The local change of the sample surface in line milling dynamically alters the milling conditions – an effect that can be exploited to speed up material removal.

While line milling is potentially faster, redeposition fills up most of the trench and hence restricts the viewable cross-section. Targeting a specific cross-section depth requires experimentation and can be cumbersome. With Fastmill, a newly-introduced scanning strategy, milling speed is enhanced by optimally exploiting the angle-dependent sputtering effect. Fastmill enables up to 40% faster milling than regular line milling. You need only activate one checkbox in the regular cross-sectioning or TEM prep workflow wizards.



Comparison of milling strategies in silicon. Material removal with conventional milling (left) takes 10 min 54 sec whereas with Fastmill the same amount of material is removed in 7 min 21 sec (right).

Expand Your Possibilities


- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

ToF-SIMS enables High Throughput in 3D Analysis


Secondary Ion Mass Spectroscopy (SIMS) is an established means of analyzing surfaces that gives you excellent sensitivity and mass resolution, along with the ability to differentiate between isotopes. Adding ToF-SIMS (time of flight secondary ion mass spectrometry) to your Crossbeam brings unique analytical capabilities to the FIB-SEM.

You will Benefit from:


- parallel detection of atomic and molecular ions down to the ppm level
- analysis of light elements, e.g. lithium
- analysis of isotopes
- analytical mapping and depth profiling
- better than 35 nm lateral resolution, 20 nm depth resolution
- post-mortem retrieval of any signal from the ROI




Working principle of SIMS: The Ga focused ion beam (blue) removes material from the top few nm of the sample surface. Different sputtered ion species (light and dark grey) are collected and transferred to the ToF-SIMS detector.



Top: Al (27 amu) map of a calibrated BAM L200 sample. The FOV is 2 μm . Bottom: Line profile for the area within the green frame. Lines with a width and separation of 33.75 nm can be resolved clearly (arrows).



SIMS spectrum of boron doped silicon. The peaks at 10 and 11 amu correspond to the two isotopes of boron. The concentration of ^{10}B is below 4.2 ppm.



Left: SEM image of the cross section of an AlAs GaAs multilayer system. The AlAs layers are 10 nm thick. Right: Corresponding SIMS depth profile showing the aluminum signal at 27 amu for the top 11 layers.

Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Gemini's Novel Optics – Profit from Surface Sensitive Imaging



High resolution imaging at low landing energy has become standard in SEM applications. It is not only required for beam sensitive samples and for non-conductive materials, it is also applied for gaining true surface information without undesirable background signal from deeper sample layers.

SEM imaging performance of the Gemini optics is optimized dramatically at low and very low voltages through the introduction of a novel optical design. It includes a high resolution gun mode which results in the reduction of the primary beam energy width by 30%. And that is finally responsible for the advances in resolution.


Additionally, a two-step deceleration modus, the so-called Tandem decel, is introduced with the novel optical design of Gemini columns. The electron optical column of ZEISS FE-SEMs have an integrated beam deceleration by design using the beam booster technology. Now, an additional external sample biasing further improves low voltage resolution and contrast. A high negative bias voltage is applied to the sample, which decelerates the electrons of the primary electron beam, thus effectively reducing landing energy:

$$E_{\text{landing}} = E_{\text{primary}} - E_{\text{bias}}$$

The Tandem decel mode can be used in two different application modes: one for contrast enhancement by applying a variable negative bias voltage between 50 V and 100 V and the second enables low voltage resolution improvement by applying a negative bias voltage of 4 different fixed values of 1 kV, 2 kV, 3 kV or 5 kV.



Brain tissue sample, showing numerous nerves that are surrounded by layers of special molecules for insulation, the myelin sheaths. Imaged at 1 kV without (left) and with Tandem decel (right). With the bias activated the myelin sheaths are clearly visible. Sample courtesy of: M. Cantoni, EPFL Lausanne, CH.



Sample biasing applies a voltage of up to 5 kV, using the optional feature Tandem decel, and improves imaging with the Gemini lens at low voltages even further.



Tandem decel sample holder for 9 single specimens.

Expand Your Possibilities


- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Make the Most of Your Crossbeam


Start smart, without that time-consuming search for the region of interest on your sample: take advantage of the optional navigation camera on the airlock. Locate specimens or specific sites, even in color. The integrated user interface makes it easy to navigate to your ROI. Select the large airlock and handle wafers of up to 8-inch diameter with fast sample transfer times.

Because it can be configured with two chamber sizes, Crossbeam 550 guarantees a high level of flexibility. The large chamber lets you customize your Crossbeam with a wider range of imaging, analytical and sample modification capabilities. Opt for a multi-channel gas injection system (GIS) to inject up to five different gases or configure your Crossbeam with up to two single GIS systems.

The large chamber offers you the possibility of configuring three pneumatically-driven accessories simultaneously, e.g. an aSTEM (annular scanning transmission electron microscopy) detector, an annular backscatter detector and a local charge compensation.



The navigation camera on the airlock helps you find your region of interest quickly and easily.



Crossbeam equipped with two Uni-GIS units, both configured for optimal access angles to achieve optimal depositions.



[▶ Click here to view this video](#)

Take high resolution images in transmission mode with the STEM detector and exploit all contrast mechanisms from brightfield to high angular annular darkfield.

Tailored Precisely to Your Applications

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service

Typical Applications, Typical Samples	Task	ZEISS Crossbeam Offers
Cross-sectioning	Acquire high resolution images of cross-sections to obtain sub-surface information.	Crossbeam offers a wide range of detectors for a comprehensive characterization of your sample. Up to four detector signals can be acquired simultaneously to get more information at the same time. The Gemini lens design does not expose your sample to a magnetic field. It allows distortion free imaging of large fields of view. Coupled with image frame store resolutions of up to 50 k x 40 k pixel your Crossbeam is ideal for large area mapping applications.
FIB-SEM Tomography	Perform serial cross-sectioning to image and reconstruct volumes of your sample	The Inlens EsB detector provides excellent material contrasts and allows surface sensitive imaging because it reduces the information depth to only a few nanometers. When used during milling with the focused ion beam, it speeds up long-term experiments. Intelligent software solutions enable long and unattended tomography runs for reliable and precise results in the shortest time.
FIB-SEM Tomography in Life Sciences	Acquire high resolution images of your cross-section or perform large volume tomography for morphological analysis.	Precisely target, image and reconstruct the volume of interest to get 3D information from your biological samples.
3D Analytics	Study the chemical and crystallographic microstructure of your sample.	Crossbeam is the perfect tool for 3D EDS and 3D EBSD analysis of your sample. Different packages are provided for fully automated acquisition of the 3D datasets.
Conventional TEM Sample Preparation	Prepare thin lamellae for analysis in TEM or STEM.	Crossbeam offers a complete solution for preparing TEM lamellae, even for batches. Benefit from the low voltage performance of the Ion-sculptor FIB column in gaining high quality lamellae and avoid amorphization of delicate specimens. For conventional TEM sample preparation, use a guided, semi-automated workflow with three simple steps to get started and wait for automatic execution. For preparing high quality lamellae, use the patented X ² sample holder during final thinning. Benefit from an endpoint detection software that gives you accurate information about the thickness of your lamella.
Fully Automated TEM Sample Preparation	Achieve best sample quality and maximize lab productivity with unattended, automated TEM lamella preparation for silicon devices in the semiconductor lab environment and benefit from high reliability and success rates.	For efficiency and higher productivity with high throughput and ease of use, Crossbeam 550 Samplefab provides a unique automated platform. It allows the full lamella preparation workflow to run unsupervised for multiple TEM samples. Run automated chunking, lift-out, attach to grid, and thinning steps according to operator-defined recipes. Experience the automation yield promise of >90% automation success without user intervention from chunking to grid attach and retain flexibility for manual preparation when desired.
Nanopatterning	Create new structures or modify existing structures by FIB (or SEM) and different gases.	Perform FIB patterning tasks with full SEM control in real time. Just choose and draw the shapes you want to create on your FIB image, set up the parameters and start patterning. The system's user-friendly software helps inexperienced users to achieve great results. For most advanced fabrication tasks, the software allows you to access all relevant SEM, FIB or GIS parameters to tailor the best FIB patterning strategy at single object level. You can plan and create your FIB exposure work offline.
Surface sensitive analysis of batteries or polymers	Characterize the composition of the first few atom layers of solid surfaces.	Adding the ToF-SIMS spectrometer lets you analyze trace elements such as lithium, detect isotopes, and perform elemental mapping and depth profiling with a lateral analytical resolution down to 35 nm.
Gain speed and quality in in situ studies	Localize and access deeply buried structures quickly	Add a femtosecond laser to your Crossbeam, find hidden ROIs using the correlative workspace GUI in ZEISS Atlas 5 e.g. with XRM datasets, ablate materials quickly and safely with the femtosecond laser.

ZEISS Crossbeam at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service


Nanopatterning



Live imaging while milling a spiral: SE signal (left), Inlens SE (right).






SEM Imaging



Alumina nanospheres imaged at 1 kV and FIB-SEM coincidence point with Tandem decel exemplifies high resolution, surface sensitive imaging of challenging samples.

Cross-sectioning and 3D Analysis






LiMn_2O_4 cathode material of a lithium ion battery. Close-up of cross-section shows surface information on an Inlens SE image and unique, pure materials contrast with an Inlens EsB image. The distribution of lanthanum (red) and manganese (green) is derived from an EDS map (from left to right).

ZEISS Crossbeam at Work


- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service

Nanopatterning




Nanofluidics channels fabricated by FIB in a silicon master stamp (left). Detail: meander-shaped channel (center). Inlets and outlets have a funnel shape (right). Courtesy of: I. Fernández-Cuesta, INF Hamburg, Germany.


Milling



Trenches milled in high entropy alloy, dimensions $25 \mu\text{m} \times 15 \mu\text{m}$, milling time 3 minutes for 65 nA box (right) and 11 minutes for 30 nA box (left).



Trench milled in silicon, dimensions $100 \times 30 \times 25 \mu\text{m}^3$, milling time 10 minutes using 100 nA FIB current.




Pillar for compression testing after being machined using lathe milling.


ZEISS Crossbeam at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service

TEM Sample Preparation




SEM view of a TEM lamella chunk cut through a hard disk read sensor on ceramic substrate, ready to be lifted out.




Ion-sculptor 1 kV image of a lamella prepared from a cell phone processor. Tight beam profile at low energies allows precise, artifact-free thinning of the lamella.


LaserFIB Preparation




Pillar array in silicon, laser milled in about 30 seconds, ready for fine polishing with the Gallium FIB. Note the excellent surface finish between the pillars and the high contour accuracy and surface quality achieved by laser processing only. This will enable further processing or investigations easily e.g. mechanical tests or APT (atom probe tomography) sample preparation.



Silicon in $\langle 110 \rangle$ orientation, STEM image of a FIB lamella in a TEM. $\langle 110 \rangle$ silicon dumbbells and a twin boundary are clearly visible. The TEM lamella was prepared with the Ion-sculptor FIB of ZEISS Crossbeam 550 with low kV thinning. Image Courtesy of: C. Downing, CRANN Institute, Trinity College, Dublin, Ireland. Nion UltraSTEM 200.



Chromium carbides at the grain boundary of thermally-affected X2CrNi18-10 steel: STEM BF in a FIB-SEM (left), EDS chromium map (right).




H-bar lamella preparation by fs-laser on a copper semi-circle grid. The left lamella is 400 μm wide, 215 μm deep and has a thickness of about 20 μm at the top. It was cut by the laser in 34 s. The amount of material that needs to be removed by FIB for final thinning is significantly reduced.

ZEISS Crossbeam at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service


FIB – Tomography in Life Sciences Cell Biology – HeLa Cells



▶ [Click here to view this video](#)

Investigation of different cell compartments in single cells. Individual HeLa cells were grown in culture dishes, chemically-fixed and resin-embedded in EPON. Voxel size $5 \times 5 \times 8$ nm, Inlens EsB detection, 1400 sections. 3D visualization with Dragonfly Pro, ORS. Courtesy: A. Steyer and Y. Schwab, EMBL, Heidelberg, DE.


Cell Biology – Algae



▶ [Click here to view this video](#)

3D reconstruction of a vitrified Emiliania huxleyi coccolithophore obtained from a cryo-FIB-SEM image series. The 3D reconstruction shows the immature coccolith (in yellowish), a coccolith in statu nascendi (blue) and lipid bodies (red). Courtesy: L. Bertinetti, Max-Planck Institute of Colloids and Interfaces, Potsdam, DE and A. Scheffel, Max-Planck Institute Plant Physiology, Potsdam, DE.


Neuroscience – Brain Sections



▶ [Click here to view this video](#)

Large area milling and imaging of a brain section with the 3D module of ZEISS Atlas 5. High current allows fast milling and imaging of large fields of view up to $150 \mu\text{m}$ in width. The depicted brain image has a field of view of $75 \mu\text{m}$ in width and was milled with a beam current of 20 nA. Courtesy: C. Genoud, FMI Basel, CH.


Developmental Biology – C.elegans



▶ [Click here to view this video](#)

Understanding the morphology of a whole organism in 3D with the highest resolution and reliability. The data set shows a large 3D volume of C.elegans consisting of 10.080 z-sections at $5 \times 5 \times 8$ nm pixel size. The nematode was high pressure frozen and freeze-substituted in EPON. Even the smallest structures inside the worm can be identified very easily. Courtesy: A. Steyer and Y. Schwab, EMBL Heidelberg, DE; and S. Markert and C. Stigloher, University of Wuerzburg, DE.


Microbiology – Trypanosoma



Ultrastructural investigation of the parasite Trypanosoma brucei. The cells are high pressure frozen and freeze-substituted in EPON. Acquisition of 800 z-sections which corresponds to $\sim 8 \mu\text{m}$ thickness in z; pixel size in x/y is 5 nm. Sample courtesy: S. Vaughan, Oxford Brookes University, Research Group 'Cell biology of Trypanosomes', UK.

Your Flexible Choice of Components

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service



Available Options

1. Focused Ion Beam (FIB) column
2. Electron flood gun allows ion beam preparation of non-conductive samples
3. Local charge compensation allows SEM imaging and analysis of non-conductive samples
4. Retractable ToF-SIMS spectrometer for parallel mass detection with excellent spatial resolution
5. Multichannel Gas Injection system (GIS) for up to 5 precursor materials on a single flange
6. Uni-GIS for high angle sample access,

- two systems configurable
7. Manipulators for sample handling and probing
8. Annular STEM for high resolution transmission imaging or aBSD4 detector for high efficiency and angle selective material characterization
9. Femtosecond laser for massive material ablation
10. Airlock solution (80 mm or 200 mm sample size) for fast and efficient sample transfer and fast pumping times with integrated navigation camera

Further Options

- Inlens EsB detector for high resolution without topographic artifacts and unique material contrast
- SESI detector for secondary electron and secondary ion imaging
- Atlas 5 for advanced tomography, patterning and EDS and EBSD analytics in 3D
- Plasma cleaner
- Electrostatic beam blanker for SEM
- Tandem decel for enhanced resolution and contrast at low voltages for suitable samples
- Analytic detectors: EDS, WDS, EBSD
- 34 inch, 21:9, ultra-wide screen monitor

Technical Specifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › **Technology and Details**
- › Service

	ZEISS Crossbeam 350	ZEISS Crossbeam 550	
Advantages	Maximum sample variety due to variable pressure mode, wide range of <i>in situ</i> experiments.	High throughput in analytics and imaging, high resolution under all conditions.	
SEM	Schottky Emitter	Schottky Emitter	
	1.7 nm @ 1 kV	1.4 nm @ 1 kV	
	1.5 nm @ 1 kV with Tandem decel	1.2 nm @ 1 kV with Tandem decel	
	1.2 nm @ 1 kV / DCV*	1.1 nm @ 1 kV / DCV*	
	1.9 nm @ 200 V with Tandem decel	1.6 nm @ 200 V with Tandem decel	
	0.9 nm @ 15 kV	0.7 nm @ 15 kV	
	0.7 nm @ 15 kV / DCV*	0.6 nm @ 15 kV / DCV*	
	0.7 nm @ 30 kV (STEM mode)	0.6 nm @ 30 kV (STEM mode)	
	2.3 nm @ 1 kV (WD 5 mm)	1.8 nm @ 1 kV (WD 5 mm)	
	1.7 nm @ 1 kV with Tandem decel (WD 5 mm)	1.3 nm @ 1 kV with Tandem decel (WD 5 mm)	
	1.1 nm @ 15 kV (WD 5 mm)	0.9 nm @ 15 kV (WD 5 mm)	
	2.3 nm @ 20 kV & 10 nA (WD 5 mm)	2.3 nm @ 20 kV & 10 nA (WD 5 mm)	
	Beam current: 5 pA – 100 nA	Beam current: 10 pA – 100 nA	
FIB	LMIS: Lifetime: 3000 µAh	LMIS: Lifetime: 3000 µAh	
	Resolution: 3 nm @ 30 kV (statistical method)	Resolution: 3 nm @ 30 kV (statistical method)	
	Resolution: 120 nm @ 1 kV & 10 pA (optional)	Resolution: 120 nm @ 1 kV & 10 pA	
	Beam current: 1 pA – 100 nA	Beam current: 1 pA – 100 nA	
Detectors	Inlens SE, Inlens EsB, VPSE (Variable pressure secondary electron detector), SESI (secondary electron secondary ion), aSTEM (scanning transmission electron), aBSD (backscatter detector)	Inlens SE, Inlens EsB, ETD (Everhard-Thornley detector), SESI (secondary electron, secondary ion), aSTEM (scanning transmission electron), aBSD (backscatter detector), CL (cathodoluminescence)	
Chamber Size and Ports	Standard with 18 configurable ports	Standard with 18 configurable ports	Large with 22 configurable ports
Stage	X / Y = 100 mm	X / Y = 100 mm	X / Y = 153 mm, 202 mm (optional)
	Z = 50 mm, Z' = 13 mm	Z = 50 mm, Z' = 13 mm	Z = 50 mm, Z' = 20 mm
	T = -4° to 70°, R = 360°	T = -4° to 70°, R = 360°	T = -15° to 70°, R = 360°
	Repeatability: 500 nm in x/y at 0° tilt over full travel range (optional)	Repeatability: 500 nm in x/y at 0° tilt over full travel range (optional)	
Charge Control	Flood Gun	Flood Gun	
	Local Charge Compensation	Local Charge Compensation	
	Variable Pressure	-	

* digital resolution (deconvoluted)

Technical Specifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › **Technology and Details**
- › Service

	ZEISS Crossbeam 350	ZEISS Crossbeam 550
Gases	Uni-GIS: Pt, C, SiO _x , W, H ₂ O Multi-GIS: Pt, C, W, Au, H ₂ O, SiO _x , XeF ₂	Uni-GIS: Pt, C, SiO _x , W, H ₂ O Multi-GIS: Pt, C, W, Au, H ₂ O, SiO _x , XeF ₂
Store Resolution	512 × 384 to 32 k × 24 k (up to 50 k × 40 k with optional Atlas 5 3D Tomography module)	512 × 384 to 32 k × 24 k (up to 50 k × 40 k with optional Atlas 5 3D Tomography module)
Analytic Options	EDS, EBSD, WDS, SIMS, others on request	EDS, EBSD, WDS, SIMS, others on request
Retractable ToF-SIMS Spectrometer		
Detection Limit	< 4,2 ppm boron in silicon	
Lateral Resolution	< 35 nm	
Mass/Charge Range	1-500 Th	
Mass Resolution	m/Δm > 500 FWTM	
Depth Resolution	<20 nm AlAs/GaAs multilayer system	
Femtosecond Laser		
Type	Diode pumped solid-state laser (cat.1)	
Wavelength (λ)	515 nm (green)	
Optics	telecentric	
Pulse Duration	<350 fs	
Spot Size	<15 μm	
Automated TEM Sample Preparation on Crossbeam 550 Samplefab (full specification: www.zeiss.com/samplefab)		
Advantages	Your robust solution for unsupervised, automated TEM lamella preparation	
Sample manipulation	Integrated manipulator: Oxford OmniProbe 400, 4 axes including rotation	
Workflow	Fully automated multi-site workflow: 1. Chunking 2. Lift-out (with transfer to grid) 3. Thinning to 100 nm	
Automation Yield	Maximum of 1 in 10 lamellae require operator intervention during chunking & lift-out (steps 1,2)	
Throughput	10 lamellae in 8 hours for the complete process (steps 1,2,3)	
Endpoint Control	Continuous SEM imaging, simultaneously with FIB milling, and on-the-fly fine adjustment of milling for manual endpoint control	

ZEISS Service – Your Partner at All Times

Your microscope system from ZEISS is one of your most important tools. For over 175 years, the ZEISS brand and our experience have stood for reliable equipment with a long life in the field of microscopy. You can count on superior service and support - before and after installation. Our skilled ZEISS service team makes sure that your microscope is always ready for use.

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › Technology and Details
- › **Service**

Procurement

- Lab Planning & Construction Site Management
- Site Inspection & Environmental Analysis
- GMP-Qualification IQ/OQ
- Installation & Handover
- IT Integration Support
- Startup Training

Operation

- Predictive Service Remote Monitoring
- Inspection & Preventive Maintenance
- Software Maintenance Agreements
 - Operation & Application Training
 - Expert Phone & Remote Support
- Protect Service Agreements
 - Metrological Calibration
 - Instrument Relocation
 - Consumables
 - Repairs

New Investment

- Decommissioning
- Trade In

Retrofit

- Customized Engineering
- Upgrades & Modernization
- Customized Workflows via ZEISS arivis Cloud



Please note: Availability of services depends on product line and location

www.zeiss.com/microscopy/service



Carl Zeiss Microscopy GmbH
07745 Jena, Germany
microscopy@zeiss.com
www.zeiss.com/crossbeam

Follow us on social media:

