SYCOP 3

System Control Panel for Zoom and Stereo Microscopes

**Quick Reference Guide
Kurz-Gebrauchsanweisung
Mode d'emploi succinct
Manual breve
Guia Rápido de Instruções de Uso
Guida breve**

English

Knowledge of this manual is required for the operation of the instrument. Therefore it's mandatory to make yourself familiar with the contents of this manual and pay special attention to hints concerning safe operation of the instrument. The specifications are subject to change; the manual is not covered by an update service. Unless expressly authorized, forwarding and duplication of this document and the software package, as well as utilization and communication of its contents are not permitted. Violations will entail an obligation to pay compensation. All rights reserved in the event of granting of patents or registration of a utility model. All names of companies and products mentioned in this manual may be trademarks or registered trademarks. Quoting of product names is for information only and does not represent any trademark misuse. Carl Zeiss Microscopy GmbH is not liable for the performance or the use of these products.

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435611-6044-008
Version 1, 09/30/2017

1 PRELIMINARY NOTE

This quick reference guide contains important information on device safety, the intended use of the device and how to use it safely. The complete operating manual contains further information on usage, maintenance, changing fuses and troubleshooting.

The complete SYCOP 3 operating manual is available for download at www.zeiss.com/AxioZoomV16

2 COMPLIANCE NOTICE

SYCOP 3 and EMS 3:

have been designed, produced and tested in compliance with

DIN EN 61010-1 (IEC 61010-1) "Safety requirements for electrical equipment for measurement, control and laboratory use"

and meet the requirements of the

"RoHS Directive" 2011/65/EG

also meet the

DIN EN 61010-2-101 "Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-101: Special requirements for In Vitro Diagnostic (IVD) medical devices"

as well as the

98/79/EC Appendix 1 for IVD products (In Vitro Diagnostic)

Instrument identification marking

 mark

Instrument disposal according to

EU Directive 2012/19/EC (WEEE).

3 NOTES ON INSTRUMENT SAFETY

Particular attention must be paid to the following warnings:

| | |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | <p>DANGER</p> <p>Risk of electric shock</p> <ul style="list-style-type: none"> – The power plug is equipped with an protective earth conductor, therefore, it has to be used with an appropriately equipped socket. – The protective capacity must not be rendered ineffective by the use of extension cables with no earth conductor. – In case of error there is a risk of contact with life-threatening mains voltage. |
|  | <p>DANGER</p> <p>Risk of damage to the device from penetrating moisture</p> <ul style="list-style-type: none"> – Devices with a mains connection and their accessories (particularly external table power pack) must not come into contact with moisture or liquids. – If the housing of devices or accessories that conduct mains voltage (e.g. table power pack) is damaged, these must be permanently taken out of service and disconnected from the mains. – Non-compliance in both cases results in exposure to life-threatening mains voltage. |
|  | <p>WARNING</p> <p>Risk of electric shock</p> <ul style="list-style-type: none"> – Always switch off the device and pull out the power plug before cleaning or opening the device. – Disconnect the power plug before changing device fuses. Failure to do so may result in contact with life-threatening mains voltage. |
|  | <p>CAUTION</p> <p>Risk of electric shock, overheating of the power cable, in worst case, fire</p> <p>Where there are detachable power cables, please note the following:</p> <ul style="list-style-type: none"> – Detachable power supply cables must not be replaced with power supply cables that do not meet the specifications. – Use only the specified power cables to avoid these becoming severely overheated. |

| | |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | <p>CAUTION</p> <p>Risk of ineffective safety measures</p> <ul style="list-style-type: none">– If it is determined that protective measures are no longer effective, the device must be taken out of service and secured against inadvertent operation, to avoid the risk of injury.– Contact ZEISS Service to have the device repaired. |
|  | <p>CAUTION</p> <p>Risk of fire from heat accumulation</p> <ul style="list-style-type: none">– Always check the device is switched off before covering it with a dust cover.– Clogging or covering of ventilation slats is also prohibited.– Failure to comply with these two instructions can lead to an accumulation of heat, which can damage the device and, in the worst-case scenario, cause a fire. |
|  | <p>CAUTION</p> <p>Risk of damage to the device due to incorrect or improvised fuses</p> <ul style="list-style-type: none">– When changing the device fuses, only fuses specified in the technical data may be used.– The use of other fuses and short-circuiting the fuse box are prohibited, to avoid the risk of overheating, smoke formation or fire in the device. |
|  | <p>NOTICE</p> <p>Environmental damage due to improper disposal</p> <ul style="list-style-type: none">– Do not dispose of defective equipment in regular domestic waste. It should be disposed of in accordance with the prevailing national legal requirements. The manufacturer of the device is under the legal obligation to take defective devices back.– Sample waste should be disposed of in accordance with the prevailing national legal requirements. |

4 NOTES ON WARRANTY

The manufacturer guarantees that the instrument is free of material or manufacturing defects upon delivery. Any defects must be reported immediately and steps taken to minimize damage. If such a defect is reported, the instrument manufacturer shall be obliged to correct the fault, either by repairing the instrument or replacing it with a new one, at the manufacturer's discretion. No warranty is given for defects caused by natural wear and tear (particularly of wearing parts) and improper use of the instrument.

The instrument manufacturer shall not be liable for damage caused by misuse, negligence or any other tampering with the instrument, particularly the removal or replacement of instrument components, or the use of accessories from other manufacturers. Such actions shall invalidate any warranty claims.

With the exception of the work described in this quick reference guide, no maintenance or repair work is to be carried out on the SYCOP 3. Repairs may only be performed by ZEISS Service or individuals specially authorized by ZEISS Service. In the event of a problem with the instrument, please contact your local ZEISS representative.

5 INTENDED USE

With the SYCOP 3 system control panel and the EMS 3 controller required for it, the following stereo and zoom microscopes can be operated:

- SteREO Discovery.V12/V20,
- SteREO Lumar.V12,
- Axio Zoom.V16.

The correct operation of the above microscope systems is essential for their safe and successful functioning. Before setting up, putting into operation and utilizing the SYCOP 3, be sure, therefore, to carefully read this Operating Manual and the relevant operating manuals for the microscopes connected and for other components and accessories that may be used with them and to make yourself thoroughly familiar with their contents.

Keep the operating manuals in a place where all the users can, at any time, access the information required for the use of the microscope system.

CAN components connected to the system are controlled and monitored via joystick, scroll wheels, pushbuttons, (optional) footswitches and touchscreens.

6 ADDITIONAL APPLICABLE OPERATING MANUALS

In addition to the present operating manual, the following operating manuals for optional devices should also be consulted depending on the equipment of the system:

- ZEN Software Description (blue edition) (online version)
- Operating manual for SteREO Discovery V12/V20
- Operating manual for Axio Zoom.V16
- Operating manual for the computer system
- Operating manual for the control computer
- Operating manual for the monitor
- Operating manual for the microscope stage

7 TECHNICAL DATA

Dimensions and Weight

| | |
|----------------------------------------|----------------------------------|
| SYCOP 3 (length x depth x height)..... | 160 mm x 260 mm x 125 mm; 0.6 kg |
| EMS 3 (length x depth x height)..... | 250 mm x 220 mm x 105 mm; 3 kg |

Power supply

| | |
|--------------|-----------------------------------|
| SYCOP 3..... | Via MDR 2x20 cable from EMS 3 |
| EMS 3..... | 100 V AC to 240 V AC, 50 Hz/60 Hz |

Operating Data

| | |
|--------------------------------------|---------------------------------------------------------------------------------------------|
| Operating environment | Enclosed spaces |
| Safety class | I |
| Degree of protection..... | IP 20 |
| Electrical safety | As specified in DIN EN 61010-1 (IEC 61010-1), taking into account CSA and UL regulations |
| Overvoltage category | II |
| Radio interference suppression | As specified in EN 55011 Class A |
| Noise immunity..... | As specified in DIN EN 61326-1 |
| AC line voltage range..... | 100 V AC to 240 V AC \pm 10 % |
| Line frequency | 50 Hz to 60 Hz |
| Power consumption of EMS 3 | 180 VA |
| Fuse protection of EMS 3 | 2x T 4 A/H 250 V |

Ambient Conditions

Shipping (in Packaging)

| | |
|---------------------------------------|---------------|
| Permissible ambient temperature | -40 to +70 °C |
|---------------------------------------|---------------|

Storage (in Packaging)

| | |
|----------------------------------------------------------|---------------|
| Permissible ambient temperature | -40 to +70 °C |
| Maximum permissible humidity (without condensation)..... | 75 % at 35 °C |

Operation

| | |
|----------------------------------------------------------|----------------------|
| Permissible ambient temperature | +5 °C to +35 °C |
| Maximum permissible humidity (without condensation)..... | 75 % at 35 °C |
| Operating environment | Enclosed spaces |
| Pollution degree | 2 |
| Maximum operating altitude..... | 2,000 m |
| Atmospheric pressure | 800 hPa to 1,060 hPa |

8 INSTALLATION CONDITIONS

The installation conditions for SYCOP3 correspond to those of the microscopes assigned to them: SteREO Discovery.V12/V20, SteREO Lumar.V12, Axio Zoom.V16.

The ambient conditions (humidity, temperature, pressure) specified in Section 7 should be adhered to.

When documenting fluorescence preparations, do not expose the microscope to direct sunlight. This helps minimize the influence of stray light.

For prolonged experimental procedures using the microscope (e.g. involving incubation or long image-recording procedures) choose an installation location with minimum vibration.

Vibration stemming e.g. from forced ventilation, floor impact noise or other shocks can have an impact on the supplementary modules of the microscope, some of which react sensitively. It can also lead to reduced image quality and experimental errors.

To reduce vibrations at the microscope installation location, ZEISS offers various passive and active damping options (Axio Imager (000000-0477-190 anti-vibration plate); small, air-damped system table (000000-1984-812), table top mat TS 150).

Select an installation site which is as temperature-stable as possible for long-term observations.


For reliable long-term observations, the device must first warm up. The temperature equilibrium is attained within 2-3 hours, depending on the light source, electronics and room temperature control. This time must be planned into critical studies.

Space requirements

The distance of the system to the wall should be at least 15 cm, in order to ensure sufficient air circulation and accessibility of the cabling..

9 ESSENTIAL OPERATION CONTROLS

Operation and Function Controls on the SYCOP 3 and on the EMS 3

- 1 Joystick
- 2 Scroll wheel
- 3 Scroll wheel
- 4 Pushbutton
- 5 Pushbutton
- 6 Pushbutton
- 7 Pushbutton
- 8 Pushbutton
- 9 Pushbutton
- 10 Touchscreen
- 11 Display line
- 12  Button (Standby)
- 13 Pilot lamp
- 14 **Power** switch (I/O)

Fig. 1 Operation and Function Controls on the SYCOP 3 and the EMS 3

The **SYCOP 3** system control panel combines four controls in one unit:

- Joystick (Fig. 2/1) permitting the two principal functions of the microscope to be controlled in a combined way: zooming and focusing
- Scroll wheels (Fig. 2/2) facilitating the fine focusing of the microscope image in greater magnifications
- Pushbuttons (Fig. 2/3) arranged in six pairs to allow additionally selectable microscope functions to be individually controlled
- Touchscreen (Fig. 2/4) allowing specific microscope functions to be controlled and the overall system to be monitored

The SYCOP 3 has been designed for one-hand operation of the microscope. This allows the operator to use all the principal functions without having to interrupt viewing through the eyepieces.

The SYCOP 3 can be placed to the right or left of the microscope, independent of the location.

The EMS 3 controller is required for the SYCOP 3.


Fig. 2 Overview of Operation Controls on the SYCOP 3


Fig. 3 Initialization Screen on the Touchscreen


Fig. 4 Stage Initialization Query Window

Fig. 5 'Home' Home Page (System Unconfigured)

10 SWITCHING ON SAFELY

Check that the microscope system has been correctly set up and connected.

All transport locks are properly removed.

Switch on all the light sources and illuminators existing in the system. In so doing, follow the separate operating instructions of each device.

Now, all components can be switched on according to the operating manual of the main device.

- Before switching the device on, check that the joystick on the SYCOP 3 is in zero position and that the footswitch, if any, is freely accessible.
- Set the **Power** switch at the back of the EMS 3 controller to position **I**. If necessary, briefly press the **⏻** button (standby) on the front of the EMS 3 controller.

The system will be switched on and will initialize (Fig. 3).

If a motorized stage is connected, a stage initialization query will appear (Fig. 4).

- If the **Yes** button is pressed, initialization will be performed immediately.
- If the **No** button is pressed, initialization will be stopped and may be carried out later.

When the initialization is completed, the **Home** home page (Fig. 5) will appear.

The system is ready for subsequent configuration.

Further options for setting user-dependent parameters can be found in the main operating manual.

11 OPERATING THE SYCOP 3

11.1 Joystick

In the default setting, the joystick (Fig. 2/1) is used to control the two principal functions of the microscope (zooming and focusing).

Slight, careful deflection of the joystick results in fine zooming and focusing movements, while full deflection produces faster zooming and focusing movements.

| | | |
|---------------------------------------|----------|--------------------------------------------|
| - Tilt the joystick <i>left</i> . | Zooming | Reduce zoom magnification. |
| - Tilt the joystick <i>right</i> . | Zooming | Increase zoom magnification. |
| - Tilt the joystick <i>forward</i> . | Focusing | Move the motorized focusing unit upward. |
| - Tilt the joystick <i>backward</i> . | Focusing | Move the motorized focusing unit downward. |

Whenever one of the two motorized mechanical or measuring stages intended to be used is connected to the system, the stage may additionally be moved in X or Y direction by the joystick. If you wish to switch from microscope control (zooming and focusing) to stage control and back, briefly press the joystick.

Slight, careful deflection of the joystick will result in fine stage movements, while full deflection will produce faster stage movements.

| | |
|---------------------------------------|--------------------------------|
| - Tilt the joystick <i>left</i> . | - Stage moves in -X direction. |
| - Tilt the joystick <i>right</i> . | - Stage moves in +X direction. |
| - Tilt the joystick <i>forward</i> . | - Stage moves in +Y direction. |
| - Tilt the joystick <i>backward</i> . | - Stage moves in -Y direction. |

11.2 Scroll Wheels

In the default setting, the two scroll wheels (Fig. 2/2) are used for fine focusing.

| | | |
|---------------------------------------------------------|---------------|---------------------------------------------|
| - Turn the left or right scroll wheel <i>forward</i> . | Fine focusing | Moves the motorized focusing unit upward. |
| - Turn the left or right scroll wheel <i>backward</i> . | Fine focusing | Moves the motorized focusing unit downward. |

Whenever a motorized mechanical or measuring stage is connected additionally to the system, the scroll wheels permit additional fine movement of the stage in X or Y direction. If you wish to switch from fine focusing to stage control and back, briefly press one of the two scroll wheels.

| | |
|-------------------------------------------------|--------------------------------------------|
| - Turn the left scroll wheel <i>forward</i> . | Fine movement of the stage in -X direction |
| - Turn the left scroll wheel <i>backward</i> . | Fine movement of the stage in +X direction |
| - Turn the right scroll wheel <i>forward</i> . | Fine movement of the stage in +Y direction |
| - Turn the right scroll wheel <i>backward</i> . | Fine movement of the stage in -Y direction |

11.3 Pushbuttons

Further different microscope functions can be controlled with the three pairs of pushbuttons (Fig. 2/3). These need to be configured accordingly beforehand (see main operating manual).


11.4 Touchscreen

Aside from controlling selected microscope functions, the touchscreen (Fig. 2/4) also serves to display all the microscope settings and parameters required.

The display line on the lower edge of the touchscreen shows the current assignment of such control elements as the joystick, the scroll wheels or the pushbuttons.

If this display line is briefly touched, the following will be shown or hidden:

- the functions of the pushbuttons
- the functions of the joystick and scroll wheels.

12 SWITCHING OFF SAFELY

- Press the **Shut Down** button in the **Home** home page on the touchscreen.

When the safety query is answered by pressing the **Yes** button (Fig. 6), the system, including light sources and illuminators, will be shut down (Fig. 7). It will then be in the standby mode.

- Set the **Power** switch at the back of the EMS 3 controller to position **O**. The system will be switched off completely.
- Then switch off light sources and illuminators separately.


Fig. 6 Shut Down Safety Query Window

Fig. 7 Shutdown Window

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