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Knowledge of this manual is required for the operation of the instrument. Would you therefore please make yourself familiar with the contents of this manual and pay special attention to hints concerning safe operation of the instrument.

The specifications are subject to change; the manual is not covered by an update service.

© Unless expressly authorized, forwarding and duplication of this document, as well as utilization and communication of its contents are not permitted. Violations will entail an obligation to pay compensation.

All rights reserved in the event of granting of patents or registration of a utility model.

Herausgeber / Issued by: Carl Zeiss Microscopy GmbH
Carl-Zeiss-Promenade 10
07745 Jena, Germany

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Safety Instructions

The following warning and information symbols are used in this manual:

CAUTION

This symbol indicates a possible hazard to the user of the instrument.

ATTENTION

This symbol indicates a possible hazard to the device or system.

NOTE

This symbol designates a note that should be closely followed.

The use of the Vertical Illumination Devices S is permitted only for the purposes described below. Strictly observe also the operating manual of the microscope used (SteREO Discovery).

Dirt and dust may impair the performance of the device. It must therefore be protected from such influences to the greatest possible extent.

Do not operate the devices in potentially explosive areas. Use them in closed rooms only.

Defective Vertical Illumination Devices S must not be disposed of with the household waste; dispose of them in compliance with the provisions of the law.

Also specimens have to be disposed of properly according to the applicable legal provisions and internal work instructions.

Description

The Vertical Illumination Device S is a very suitable accessory to illuminate objects extended in depth, such as bores, recesses and holes, due to the small angle between illumination and observation and the concentration of the light towards the field center.

The Vertical Illumination Device S can be operated in combination with an existing standard cold-light source, preferably with the CL 1500 ECO. An optional, flexible light guide of 4.5 mm diameter is used to connect the Vertical Illumination Device S with the cold-light source.

It can be combined with the objectives Plan S 1,0x, Achromat S 1,0x, Achromat S 0,63x, Achromat S 0,5x, PlanApo S 1,0x and PlanApo S 0,63x.

Technical Data

Dimensions (diameter x height) 79 x 24.5 mm
Weight 0.120 kg

Ambient Conditions

Transportation (in packaging):

Admissible ambient temperature -40 to +70 °C

Storage:

Admissible ambient temperature +10 to +40 °C

Admissible air humidity (no condensation) max. 75 % at 35 °C

Operation:

Admissible ambient temperature +10 to +40 °C

Admissible relative humidity (no condensation) max. 75 % at 35 °C

Operation altitude max. 2000 m


Atmospheric pressure 800 hPa to 1060 hPa

Pollution degree 2

Environment for use closed rooms

Putting into Operation

- Attach the Vertical Illumination Device S (Fig. 1/2) from below to the objective (Fig. 1/1) of the stereomicroscope, orientate it to north-south (push-pull rod (Fig. 1/5) projects forward, i.e. to the user) and fix it using both screws (Fig. 1/6).
- Introduce the flexible light guide (Fig. 1/3) into the mount of the Vertical Illumination Device S as far as it will go (if necessary, loosen the knurled screw (Fig. 1/4) slightly) and tighten the knurled screw (Fig. 1/4).
- Connect the flexible light guide with the cold-light source.
- Switch on the cold-light source.

- 1 Objective
- 2 Vertical Illumination Device S
- 3 Flexible light guide
- 4 Knurled screw
- 5 Push-pull rod
- 6 Screw (2x)


Fig. 1 Attaching the Vertical Illumination Device S

Operation

The illumination mirror of the Vertical Illumination Device S is movable between 0 and 19 mm in north-south direction and can be locked in any position within this range. With a Plan S 1,0x objective, this corresponds to an angle of 0 to 8°.

This allows the position of the illumination mirror to be optimized in dependence on the objective and zoom setting in order to obtain the appropriate illumination of the specimen and exclude vignetting.

- For an optimum observation of the object, adjust the illumination angle of the mirror (Fig. 2/1) by shifting the push-pull rod (Fig. 2/3).
- Fix the position of the push-pull rod by tightening the clamping knob (Fig. 2/2).
- If the objective is attached to the objective slider S/doc or to the triple objective nosepiece S/doc, 6x cod. in mono position, the push-pull rod must be re-adjusted to optimize the observation of the object. If the mono position on the objective nosepiece S/doc is used, also the Vertical Illumination Device S must be re-adjusted in north-south direction.

- 1 Mirror
- 2 Clamping knob
- 3 Push-pull rod with scale

Fig. 2 Operating the Vertical Illumination Device S

The scale with 2 mm intervals on the push-pull rod allows the reproducible adjustment of the mirror position.