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Knowledge of this manual is required for the operation of the instrument. Would you therefore please make yourself familiar with the contents of this manual and pay special attention to hints concerning safe operation of the instrument.

The specifications are subject to change; the manual is not covered by an update service.

- © Unless expressly authorized, forwarding and duplication of this document, as well as utilization and communication of its contents are not permitted. Violations will entail an obligation to pay compensation.

All rights reserved in the event of granting of patents or registration of a utility model.

Herausgeber / Issued by: Carl Zeiss Microscopy GmbH
Carl-Zeiss-Promenade 10
07745 Jena, Germany

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1 INTRODUCTION

1.1 Safety Instructions

The VIS-LED has been designed, produced and tested in compliance with the DIN EN 61010-1 (IEC 61010-1) standard "Safety Requirements for Electrical Measuring, Control and Laboratory Equipment".

The device complies with the EC Directives 2006/95/EC "Low Voltage Directive" and 2004/108/EC "Electromagnetic Compatibility".

The device carries the  mark.

The device has to be disposed of in compliance with the WEEE Directive 2002/96/EC.

The present operating manual includes information and warnings that must be observed by the user.

CAUTION

Non-observance of these safety instructions involves a possible hazard to the user.

CAUTION

Hot surface!

CAUTION

LED risk group 1 according to IEC 62471, emission of LED radiation.
Do not look into the LED beam. Eyes may be damaged.

Never look directly into the light beam - neither with nor without optical instruments, even not if you want to observe the specimen. Eyes may be damaged in case of non-observance!

Warning of possible laser radiation!

Only Carl Zeiss Service personnel is authorized to remove from the VIS-LED-Laser the laser radiation warning labels and the safety screws located beneath them (no customer interface)! If this precaution is not observed, eyes may be damaged when coupling the laser into the VIS-LED.

Warning of possible laser radiation!

Only Carl Zeiss Service personnel is authorized to loosen the safety screw on the microscope (no customer interface)! If this precaution is not observed, eyes may be damaged when coupling the laser into the microscope and removing the VIS-LED.

There is burn danger at the VIS-LED lamp housing during and up to 15 minutes after operation.

Operation of the device in potentially explosive areas is not permitted.

Strictly observe also the operating instructions of the microscopes and laser systems used!

1.2 Warning Labels and Service Interfaces

VIS-LED (423053-9030-000) for classical light microscopy


Fig. 1 Warning labels on the VIS-LED (423053-9030-000)

VIS-LED-Laser (423053-9060-000) for laser systems



Fig. 2 Warning labels on the VIS-LED-Laser (423053-9060-000)

CAUTION

Only Carl Zeiss Service personnel is authorized to remove from the VIS-LED-Laser the laser radiation warning labels and the safety screws located beneath them (no customer interface)! If this precaution is not observed, eyes may be damaged when coupling the laser into the VIS-LED.

Coupling point between laser microscope and VIS-LED

1 Safety screw

Fig. 3 Coupling point between laser microscope and VIS-LED

**CAUTION**

Only Carl Zeiss Service personnel is authorized to loosen the safety screw (Fig. 3/1) on the microscope (no customer interface)! If this precaution is not observed, eyes may be damaged when coupling the laser into the microscope and removing the VIS-LED.

1.3 Disposal

Please observe the following safety notes for the disposal of the VIS-LED:

Defective devices must not be disposed of with the household waste; dispose of them in compliance with the provisions of the law.

The manufacturer of the device is under the legal obligation to take back defective devices.

1.4 Warranty Notes

The VIS-LED must not be used for techniques other than those described in these operating instructions. The manufacturer cannot assume any liability for other applications.

The manufacturer guarantees that the device is free from material or manufacturing defects when delivered. Possible defects must be notified to us immediately and steps be taken to minimize damage. If notified of such a defect, the manufacturer of the device is obligated to rectify it at its discretion, either by repairing the device or delivering an intact replacement. No guarantee is provided for defects caused by natural wear (wearing parts in particular) and improper use.

The manufacturer of the device shall not be liable for damage caused by faulty operation, negligence or any tampering with the device, particularly the removal or replacement of device components, or the use of accessories from other manufacturers. This shall lead to a forfeit of all warranty claims.

No maintenance or repair work other than the operations mentioned in the present instructions may be carried out to the VIS-LED. Repairs are only permitted to Carl Zeiss Service personnel or to persons specifically authorized by it. In case of fault please turn first to the Carl Zeiss representative responsible for your region.

2 DESCRIPTION OF THE DEVICE

2.1 Name and Intended Use

Product name: **VIS-LED**
VIS-LED-Laser

The VIS-LED is an illumination device for microscopes. It can be used for both transmitted-light and reflected-light illumination.

2.2 Technical Data

Dimensions (diameter x height) 70 x 128 mm
Weight..... 1.08 kg

Environment for useclosed rooms
Degree of protectionIP 20
Electrical safety.....in compliance with DIN EN 61010-1 (IEC 61010-1)
including CSA and UL directives
Radio interference suppression..... in compliance with EN 55011 class B
Noise immunityin compliance with DIN EN 61326

LED risk group 1 in compliance with IEC 62471
Wavelength.....400 to 700 nm, peak at 460 nm

Power supply..... through microscope or external power supply unit
Voltage 0 to 12 V DC
Max. power consumption 20 W

3 INSTALLATION AND OPERATION**CAUTION**

Only Carl Zeiss Service personnel is authorized to attach and detach the VIS-LED-Laser (423053-9060-000) for laser systems (no customer interface)! If this precaution is not observed, eyes may be damaged when coupling the laser into the microscope and removing the VIS-LED.


Fig. 4 Attaching the VIS-LED to the Axio Imager transmitted-light mount

3.1 Attaching the VIS-LED (423053-9030-000)**CAUTION**

Avoid looking directly into the LED light.

**CAUTION**

Attach and detach the VIS-LED to and from the transmitted-light mount only when it is switched off.

- Loosen the clamping screw (Fig. 4/1) on the transmitted-light mount (Fig. 4/2) or the reflected-light mount and remove the halogen lamp.
- Push the VIS-LED (Fig. 4/3) with its dovetail into the mount (Fig. 4/2) and fix it with the clamping screw (Fig. 4/1) using the ball-headed screwdriver SW 3.
- Connect the 3-pole lamp plug (Fig. 4/5) to the lower 3-pole 12 V/100 W port (Fig. 4/6) for transmitted light (TL) or the upper port for reflected light on the rear of the device or on an external power supply unit.
- If you use a coded Axio Imager stand, turn the reflected-/transmitted-light change-over switch to the transmitted-light position (TL). If you use a motorized stand, the reflected-/transmitted-light change-over is carried out over the touchscreen on the TFT display.

If the AxioCam MRm is used for documentation, the 12 % neutral-density filter included in the delivery must be used to reduce the illumination intensity (see section 3.3).

The VIS-LED is equipped with a BNC connection that allows control via TTL trigger.

3.2 Adjusting the VIS-LED

- Switch on the VIS-LED through the microscope or an external power supply unit.
- Swing the transmitted-light diffusing plate out of position and remove the transmitted-light filter wheels.
- Focus onto the specimen using an objective $\leq 40\times$ and find a free specimen area.
- Remove the eyepiece and center the illumination image within the pupil image using the adjusting screws (Fig. 5/1 and 2).
- Optimize homogeneous illumination of the pupil image using the focusing screw (Fig. 5/3).
- Swing the transmitted-light diffusing plate again into position and insert the transmitted-light filter wheels.


Fig. 5 Adjusting the VIS-LED

3.3 Inserting the Neutral-Density Filter into the VIS-LED

CAUTION

Only Carl Zeiss Service personnel is authorized to insert and remove the neutral-density filter into and from the VIS-LED-Laser (423053-9060-000) for laser systems (no customer interface)! If this precaution is not observed, eyes may be damaged when coupling the laser into the microscope and removing the VIS-LED).

If the AxioCam MRm is used for documentation, the 12 % neutral-density filter included in the delivery must be inserted to reduce the illumination intensity. Please proceed as follows:

- Detach the VIS-LED from the microscope (see section 3.1).
- Loosen the two screws (Fig. 6/1) on the VIS-LED housing.
- Remove the lamp module from the collector.
- Remove the O-ring (Fig. 6/3) located in front of the diffusing glass.
- Insert the 12 % neutral-density filter (Fig. 6/2) in front of the diffusing glass.
- Insert the O-ring (Fig. 6/3).
- Re-assemble the lamp module and the collector.
- Tighten the two screws on the VIS-LED housing.


Fig. 6 Inserting the neutral-density filter

4 MAINTENANCE, TROUBLE-SHOOTING AND SERVICE

4.1 Trouble-Shooting

Problem	Cause	Trouble-Shooting
VIS-LED not lighting	Lamp plug (Fig. 4/5) not inserted.	Insert lamp plug (Fig. 4/5).
	Microscope or external power supply unit does not supply voltage.	Check the microscope power supply or the external power supply unit. Switch on the device.
	VIS-LED may be defective.	Have the VIS-LED checked and, if necessary, replaced by the Technical Service.
	VIS-LED is not configured (for motorized stands only).	Configure the VIS-LED via the TFT-display
VIS-LED is flickering, light intensity is not stable.	Power cable not properly installed.	Connect the power cable properly.

4.2 Maintenance and Service

The VIS-LED has no parts to be serviced by the customer.

Any intervention in mechanical, optical or electronic components in the interior of the device as well as any work on its electrical system may only be performed by Carl Zeiss Service personnel or specially **authorized** personnel.

If service is required, please contact the regional representative responsible for your area, or

Carl Zeiss Microscopy GmbH
Carl-Zeiss-Promenade 10
07745 Jena, Germany

microscopy@zeiss.com
www.zeiss.com/microscopy

Carl Zeiss Microscopy GmbH
Königsallee 9-21
37081 Göttingen, Germany