ZEISS Predictive Service

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Connectivity Summary

- Windows service
ZeissPredictiveServiceAgent.exe
- Only outbound https, 256-bit TLS/SSL secured connection on port 443
- No open ports for inbound connection
- Proxy server support
- Typically no changes to firewall settings
- Microsoft Azure cloud infrastructure certified and audited according to international standards (ISO 27001, HIPAA, FedRAMP, SOC1 and SOC2)

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Condition Monitoring Summary

- Only log file data from specific log file storage locations are transmitted (see table 1). These files are fully transparent and plain-text readable
- No user specific information or files are included in continuous monitoring

Device type	Log file name	Description	Path
Electron microscopes	EM Server.log	Information about EM Server	C:\ProgramData\Carl Zeiss\SmartSEM\Log\EM Server.log
Electron microscopes	Log Folder	Includes SmartSEM Logs (GIS, FIB, ...)	C:\ProgramData\Carl Zeiss\SmartSEM\Log
Light microscopes	remote service.log.xml	Already parsed information from ZEN.log.xml and ZENCore.log.xml	C:\ProgramData\Carl Zeiss\Remote Service\remote service.log.xml
Light microscopes	ZEN.log.xml	Information about ZEN software	C:\ProgramData\Carl Zeiss\Logging
Light microscopes	ZENCore.log.xml	Information about ZEN files Core software	C:\ProgramData\Carl Zeiss\Logging
X-ray microscopes	XradiaAppn.log where n=1,2,...N used in system software v11.1 XradiaApp_n.log where n=1,2,...N used in system software v12	Information about Versa system parameters	C:\Program Files\Carl Zeiss X-ray Microscopy\Xradia Versa\XX.X.XXXX.XXXXX\Log Files
X-ray microscopes	Centre.log Comms.log FilStore.log Getter.log Healthcheck.log Serial.log TCC.log Tube.log Vac.log Warm.log	X-ray source logs	C:\Program Files\Carl Zeiss X-ray Microscopy\Xradia Versa\XX.X.XXXX.XXXXX\logs

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Technical facts

Ports	443
Protocols	HTTPS, Secure Websocket
URL of ZEISS server	predictive-service.zeiss.com www.predictive-service.zeiss.com
IP addresses	52.174.243.245
Certification of Microsoft Azure cloud infrastructure	ISO 27001, HIPAA, FedRAMP, SOC1 and SOC2
Connection security	TLS 1.2 with AES 256bit
Supported proxy settings for installation	- No authentication - Basic authentication - Digest authentication - NTLM authentication
Available and patched local area network connected to ZEISS system PC	RJ 45 connector for LAN
Security patched and measures according to your local IT security standards	e.g. installing antivirus software; please see recommendations for data security policy on ZEISS Microscopy Systems (see references)
Established and working outbound internet connection	

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