


3D structural analysis with your widefield microscope



ZEISS Apotome 3

Reliable Optical Sectioning in Fluorescence Microscopy
with Structured Illumination

zeiss.com/apotome






Seeing beyond

Reliable Optical Sectioning in Fluorescence Microscopy with Structured Illumination

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › Technology and Details
- › Service

Fluorescence microscopy is one of the most important imaging techniques in the life sciences, with high-contrast images achieved from fluorescent signals emitted by the sample. However, the objective lens also collects light from outside the focal plane. Out-of-focus light which might originate from structures above or below the focal plane, must be removed to extract the in-focus information only.

Optical sectioning with ZEISS Apotome 3 allows you to efficiently minimize out-of-focus light. It's a totally reliable way to create crisp images and 3D renderings, even of thicker specimen. Yet your widefield microscope remains just as easy to operate as always. Add Apotome Plus to get even higher contrast and distinguish structures down to 180 nm – simply brilliant optical sections.




35 µm sagittal section of adult mouse brain, imaged with ZEISS Axio Observer and ZEISS Apotome, processed with Apotome Plus. Sample courtesy of University of California, Davis/NIH NeuroMab Facility.

Simpler. More Intelligent. More Integrated.

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Brilliant optical sections under varying experimental conditions


To image structures of sizes from hundreds of micrometers to the nanometer range, you typically use objectives with different magnifications. ZEISS Apotome 3 comes with three grids of different frequencies, giving you the best resolution for each objective. You can fully focus on your experiment as ZEISS Apotome 3 adapts to your fluorophores and light source by automatically selecting the ideal grid. Enjoy high-contrast optical sections with significantly improved axial resolution compared to conventional fluorescence microscopy.



Cos7 cells (nuclei stained with Hoechst, tubulin with Alexa 488 and Phalloidin with Alexa 568) imaged with Plan Aplanachromat 63x/1.4.

Intuitive access to efficient workflows and smart data handling


With just one click in the ZEN imaging software, the advantages of optical sectioning are available for your usual image acquisition. Boost the productivity of your workflows with intelligent functions like Direct Processing: combine several processing functions and process your data while acquiring the next image. The built-in lossless data compression automatically reduces the file size, thus saving storage space and making data handling easier. Whether you want to image a single plane or conduct a complex experiment, ZEISS Apotome 3 is your choice for easy and efficient optical sectioning.



ZEN user interface: Adding Apotome Plus to a Direct Processing workflow

Structural information with resolution down to 180 nm

Compared to conventional widefield fluorescence microscopy, ZEISS Apotome 3 already significantly increases the axial resolution: you obtain brilliant optical sections that enable 3D rendering even of thick samples. With Apotome Plus, you can go one step further and achieve confocal-like image quality with your widefield microscope. The combination of structured illumination with state-of-the-art image processing substantially improves the signal-to-noise ratio. With Apotome Plus you can resolve details that were not visible before with your widefield microscope.



XZ-view of postnatal day (P) 16 mouse retina (blood vessel marker isolectin B4: Alexa Fluor 633, red; pericytic proteoglycan NG2: Alexa546, green; VE-cadherin, Alexa 488, blue). Sample courtesy of E. Montañez Miralles, University of Barcelona, Spain.


Your Insight into the Technology Behind It

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Reliable optical sections through structured illumination

ZEISS Apotome 3 positions an optical grid in the light path to create a pattern of intensity differences on your sample. If out-of-focus light is present at a certain region of the sample, the grid becomes invisible. After the fluorescence of a grid position is acquired, the grid moves to the next position. Each optical section is calculated from a number of images acquired with different grid positions. This proven, hardware-based approach built on well documented algorithms gives you the reliability that is expected from scientific data.

Apotome can be used with multi-LED and white light sources. It recognizes the magnification and automatically moves the appropriate grid into the beam path to ensure optimal imaging parameters. No matter if you want to image specimen on slides, petri dishes, or multi-well plates, Apotome generates crisp optical sections with high resolution.



Schematic illustration of the grid projection. A: Widefield image. B – D: raw images with different positions of the grid. E: resulting optical section through the sample. Out-of-focus light is efficiently removed by the structured illumination.


Your Insight into the Technology Behind It

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Image data you can trust

Compared to purely software-based approaches, Apotome uses structured illumination which has several advantages. Purely software-based methods require either prior knowledge of the sample (AI based methods) or rely on complex algorithms that have not been peer-reviewed. As software solutions can only use the acquired widefield image, users must trust that these black-box solutions exclusively produce structures that are real and do not remove structures when “enhancing” the image.

The figure shows the comparison of a widefield image, a background-subtracted image processed using a software algorithm, and an image acquired with ZEISS Apotome 3. Even though the background-subtracted image shows a high contrast that is pleasing to the eye, it falsifies information. Features are missing and structures look entirely different. Without knowing the true image, it is almost impossible to realize this. Apotome uses the information from the structured illumination combined with documented algorithms to create a crisp optical section you can trust.



Giant live fluke stained with Hoechst 33342. The homogeneous fluorescence in the inner parts of the widefield image (A) poses a serious problem for the background correction algorithms (B). Some structures remain, but generally, there are too many black spaces between the cells. This becomes visible when comparing the results to an optical section acquired with ZEISS Apotome 3 (C). Notably, the prominent rim around the structure, as seen in the background-corrected image (B), is an artifact of interference in the widefield image, which is not seen with an optical sectioning system.


Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Apotome Plus: Resolve structural information down to 180 nm with your widefield microscope.

Compared to conventional widefield microscopy, imaging with structured illumination enables more information to be obtained for the creation of optical sections.

To achieve even higher resolution and better signal-to-noise ratios, Apotome Plus combines information from your Apotome hardware with advanced image processing. Compared to purely software-based techniques, Apotome Plus is less prone to artifacts because of the additional information of the structured illumination. Use your widefield microscope to resolve structures down to 180 nm laterally without sacrificing the ease of use of Apotome. With smart Direct Processing workflows, you can process your data while image acquisition is still running. What's more, you can take the opportunity to discover previously unresolved details of your samples by applying the Apotome Plus processing to existing z-stacks acquired with Apotome 3 and ApoTome.2 systems.



Maximum intensity projection of PML nuclear bodies of U-2 OS cells transiently transfected with mRFP-SUMO-2 (magenta) and pEGFP-PML-1 (green). Note that overexpression of PML and SUMO-2 leads to an increased size of the nuclear bodies. Image sizes: $2 \times 2 \mu\text{m}$. Sample courtesy: P. Hemmerich, Leibniz Institute on Aging – Fritz Lipmann Institute e. V. (FLI), Germany.



Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service



Apotome Plus: best image quality from challenging samples

When imaging biological samples, resolution is not always most important to answer scientific questions. Aspects like signal-to-background ratio and removal of out-of-focus information are as crucial because they directly affect image quality and subsequent data analysis.

With Apotome Plus you can significantly improve the overall image quality even for dim samples with low signals. Image large volumes of tissue, cells, plants and other samples with confocal-like image quality using a camera-based microscope. No matter if you want to create well-resolved three-dimensional renderings, image thick samples with high contrast or are interested in smallest details, Apotome Plus is your choice for achieving best image quality with your widefield microscope.



Mouse spinal cord section having nucleoli stained with DAPI (blue), myelin by Olig2 (oligodendrocyte lineage cell nuclear marker) with Alexa 555 (green) and PDGFR α (oligodendrocyte progenitor cell membrane marker) with Alexa 647 (magenta). Sample courtesy of S. C. Fernández, Instituto de Investigación Sanitaria y Biomédica de Alicante (ISABIAL), Spain.





Low Signal-to-noise image of a marine vertebrate head (zoom-in) having nuclei stained with HOECST, neurons with Alexa 488 and the glia stained with CY3 and Alexa 647. Sample courtesy of A. Fuentes, Instituto de Investigación Sanitaria y Biomédica de Alicante (ISABIAL), Spain.

Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service


Apotome 3	Apotome 3 with Apotome Plus
	Creates optical sections by structured illumination
	Established, reliable technology
	Activation by moving the slider in imaging position
	Supports white light and single LED light sources
	Integrated chromatic correction
	Automatic selection of optimal grid for each objective
	Supports a broad range of magnifications and numerical apertures
	Easy to use
	Full integration in ZEN
	Automatic, fast lossless compression
	GPU data processing
	Works with ZEISS Axio Imager, ZEISS Axio Observer and ZEISS Axio Zoom.V16
Diffraction-limited resolution	Distance discrimination down to 180 nm (xy)
Confocal-like optical sectioning	Optical sectioning strength down to 460nm (50% contrast)
Good signal-to-noise ratio	Superior signal-to-noise ratio
Good background suppression	Best background suppression




Maximum intensity projection of bovine pulmonary arterial endothelial (BPAE) cells having nuclei stained with DAPI (blue), F-actin with Alexa 488 phalloidin (green) and mitochondria with MitoTracker Red CMXRos (red).

ZEISS Apotome 3 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service





Comparison of a widefield image of cortical neurons stained for DNA and microtubules.
Courtesy of L. Behrendt, Leibniz-Institute on Aging – Fritz-Lipmann-Institut e.V. (FLI), Germany.



3D rendering of a section of cortical neurons stained for DNA and microtubules. The enhanced resolution improves the image quality significantly. Courtesy of L. Behrendt, Leibniz-Institute on Aging – Fritz-Lipmann-Institut e.V. (FLI), Germany.

ZEISS Apotome 3 at Work







- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service



Autofluorescence of a Lotus Japonicus root infected with symbiotic bacteria stained with mCherry. Courtesy of F. A. Ditengou, University of Freiburg, Germany.

ZEISS Apotome 3 at Work


- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service



*Transgenic zebrafish larvae at 4 days post fertilization staining for: Glial fibrillary acidic protein, acetylated Tubulin, GFP and DNA. Embedded in 1.2% low melt agarose.
Courtesy of H. Reuter, Leibniz-Institute on Aging – Fritz-Lipmann-Institut e.V. (FLI), Germany.*

ZEISS Apotome 3 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service




Maximum intensity projection of bovine pulmonary arterial endothelial (BPAE) cells having nuclei stained with DAPI (blue), F-actin with Alexa 488 phalloidin (green) and mitochondria with MitoTracker Red CMXRos (red). Compared to the widefield image, Apotome removes out of focus light (see xz and yz view) creating a crisp optical section. Apotome Plus further improves image quality and resolves even finer structures.

ZEISS Apotome 3 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service

3D reconstruction based on optical sections acquired with Apotome and Apotome Plus




▶ [Click here to view this video](#)

Transgenic zebrafish larvae at 4 days post fertilization staining for: Glial fibrillary acidic protein, acetylated Tubulin, GFP and DNA. Embedded in 1.2% low melt agarose. Courtesy of H. Reuter, Leibniz-Institute on Aging – Fritz-Lipmann-Institut e.V.(FLI), Germany.

ZEISS Apotome 3 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service


3D reconstruction based on optical sections acquired with Apotome and Apotome Plus



Chromosomes of A. Thalina stained with REC8. Sample courtesy of S. Durand, MPI for Plant Breeding, Germany.

Your Flexible Choice of Components

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service



1 Microscope

- Axio Observer series (inverted research microscope)
- Axio Imager 2 series (upright research microscope)
- Axio Zoom.V16 (zoom microscope)
- Simple upgrading of existing systems

2 Objectives

Recommended objective classes with the highest level of image quality:

- C-Apochromat
- Plan-Apochromat
- EC Plan-Neofluar

3 Illumination

- Viluma 5, 7, 9 (LED)
- Excelitas Xylis (white light LED)
- HXP 120 C (metal halide)

4 Cameras

- Monochrome, low noise ZEISS AxioCam camera models
- Selected 3rd-party cameras

5 Accessories

- AI Sample Finder
- Definite Focus 3
- Z-Piezo
- Autoimmersion Module


6 Software

Recommended ZEN toolkits:

- Motorized Acquisition (z-stack, time lapse, Tiles & Position, Direct Processing)
- 3D (3D Visualization and Analysis, Advanced Processing)
- Connect (2D, 3D, third-party import)
- AI (Intellesis Segmentation, Intellesis Object Classification, Intellesis Denoising)

System Overview


- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service



Option:
System table small, air-cushioned, level unregulated
000000-2147-333


System Overview

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service



System Overview

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service



Technical Specifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › **Technology and Details**
- › Service

Grid Table: Apotome 3 generates optical sections of a defined thickness (in Rayleigh units, RU and microns, μm) depending on wavelength, microscope and objective used.

Data for the Use of Upright Microscopes, e.g. ZEISS Axio Imager

Objectives for Axio Imager	V	NA	Immersion	Grid/Section thickness @490 nm [RU/ μm]			DAPI with FS34	DAPI with FS49
				High grid	Medium grid	Low grid		
EC Plan-Neofluar	10x	0.3	Air	2.9/31.9	1.7/18.2	0.9/9.9	Yes	Yes
EC Plan-Neofluar	20x	0.5	Air	2.4/9.2	1.4/5.3	0.7/2.9	Yes	Yes
EC Plan-Neofluar	40x	0.75	Air	1.6/2.8	0.9/1.6	0.5/0.9	Yes	Yes
EC Plan-Neofluar	40x	1.3	Oil	2.5/2.2	1.4/1.2	0.8/0.7	Yes	Yes
EC Plan-Neofluar	63x	0.95	Air	1.0/1.1	0.6/0.7	0.4/0.4	Yes	No
EC Plan-Neofluar	63x	1.25	Oil	1.6/1.5	0.9/0.9	0.5/0.5	Yes	Yes
EC Plan-Neofluar	100x	1.3	Oil	1.0/0.9	0.6/0.5	0.4/0.3	Yes	Yes
LCI Plan-Neofluar	25x	0.8	Oil, water or glycerin	2.9/6.6	1.7/3.7	0.9/2.0	Yes	Yes
LCI Plan-Neofluar	63x	1.3	Water or glycerin	1.5/1.3	0.9/0.7	0.5/0.4	Yes	Yes
Plan-Apochromat	10x	0.45	Air	4.2/20.4	2.4/11.5	1.3/6.2	Yes	Yes
Plan-Apochromat	20x	0.8	Air	3.2/4.9	1.8/2.8	1.0/1.5	Yes	Yes
Plan-Apochromat	40x	0.95	Air	1.6/1.7	0.9/1.0	0.5/0.5	Yes	Yes
Plan-Apochromat	40x	1.3	Oil	2.5/2.2	1.4/1.2	0.8/0.7	Yes	Yes
Plan-Apochromat	40x	1.4	Oil	2.4/1.8	1.4/1.0	0.7/0.6	Yes	Yes
Plan-Apochromat	63x	1.4	Oil	1.6/1.2	0.9/0.7	0.5/0.4	Yes	Yes
Plan-Apochromat	100x	1.4	Oil	1.0/0.8	0.6/0.5	0.4/0.3	Yes	Yes
LD LCI Plan-Apochromat	25x	0.8	Oil, water or glycerin	2.9/6.6	1.7/3.7	0.9/2.0	Yes	Yes
C-Apochromat	10x	0.45	Water	4.2/20.4	2.4/11.5	1.3/6.2	Yes	Yes
C-Apochromat	40x	1.2	Water	2.2/2.0	1.2/1.1	0.7/0.6	Yes	Yes
C-Apochromat	63x	1.2	Water	1.4/1.3	0.8/0.7	0.5/0.4	Yes	Yes
LD C-Apochromat	40x	1.1	Water	2.2/2.3	1.2/1.3	0.7/0.7	Yes	Yes
Plan-Apochromat	63x	1.46	Oil	1.5/1.0	0.9/0.6	0.5/0.3	Yes	Yes
Plan-Fluar	100x	1.45	Oil	1.0/0.7	0.6/0.4	0.3/0.2	No	No
Plan-Apochromat	100x	1.46	Oil	1.0/0.7	0.6/0.4	0.3/0.2	Yes	No

Technical Specifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › **Technology and Details**
- › Service

Data for the Use of Inverted Microscopes, e.g. ZEISS Axio Observer

Objectives for Axio Observer	V	NA	Immersion	Grid/Section thickness @490 nm [RU/μm]			DAPI with FS34	DAPI with FS49
				High grid	Medium grid	Low grid		
EC Plan-Neofluar	10x	0.3	Air	2.9/31.5	1.7/18.5	0.9/9.8	Yes	Yes
EC Plan-Neofluar	20x	0.5	Air	2.3/9.0	1.4/5.4	0.7/2.9	Yes	Yes
EC Plan-Neofluar	40x	0.75	Air	1.6/2.7	0.9/1.6	0.5/0.9	Yes	No
EC Plan-Neofluar	40x	1.3	Oil	2.4/2.1	1.4/1.3	0.8/0.7	Yes	Yes
EC Plan-Neofluar	63x	0.95	Air	1.0/1.1	0.6/0.7	0.4/0.4	Yes	Yes
EC Plan-Neofluar	63x	1.25	Oil	1.6/1.5	0.9/0.9	0.5/0.5	Yes	No
EC Plan-Neofluar	100x	1.3	Oil	1.0/0.9	0.6/0.6	0.4/0.3	Yes	No
LCI Plan-Neofluar	25x	0.8	Oil, water or glycerin	2.9/6.5	1.7/3.8	0.9/2.0	Yes	Yes
LCI Plan-Neofluar	63x	1.3	Water or glycerin	1.5/1.3	0.9/0.8	0.5/0.4	No	No
Plan-Apochromat	10x	0.45	Air	4.2/20.2	2.4/11.7	1.3/6.1	Yes	Yes
Plan-Apochromat	20x	0.8	Air	3.1/4.8	1.8/2.8	1.0/1.5	Yes	Yes
Plan-Apochromat	40x	0.95	Air	1.6/1.7	0.9/1.0	0.5/0.5	Yes	Yes
Plan-Apochromat	40x	1.3	Oil	2.4/2.2	1.4/1.3	0.8/0.7	Yes	Yes
Plan-Apochromat	40x	1.4	Oil	2.4/1.8	1.4/1.1	0.7/0.6	Yes	Yes
Plan-Apochromat	63x	1.4	Oil	1.5/1.2	0.9/0.7	0.5/0.4	Yes	Yes
Plan-Apochromat	100x	1.4	Oil	1.0/0.8	0.6/0.5	0.4/0.3	Yes	No
LD LCI Plan-Apochromat	25x	0.8	Oil, water or glycerin	2.9/6.5	1.7/3.8	0.9/2.0	Yes	Yes
C-Apochromat	10x	0.45	Water	4.2/20.2	2.4/11.7	1.3/6.1	Yes	Yes
C-Apochromat	40x	1.2	Water	2.1/1.9	1.3/1.1	0.7/0.6	Yes	Yes
C-Apochromat	63x	1.2	Water	1.4/1.3	0.8/0.7	0.5/0.4	Yes	Yes
LD C-Apochromat	40x	1.1	Water	2.1/2.3	1.3/1.4	0.7/0.7	Yes	Yes
Plan-Apochromat	63x	1.46	Oil	1.5/1.0	0.9/0.6	0.5/0.3	Yes	Yes
Plan-Fluar	100x	1.45	Oil	1.0/0.7	0.6/0.4	0.3/0.2	No	No
Plan-Apochromat	100x	1.46	Oil	1.0/0.7	0.6/0.4	0.3/0.2	Yes	No

Technical Specifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › **Technology and Details**
- › Service

Apotome Plus

Resolution	Lateral resolution (x/y): After processing, distances down to 180 nm can be discriminated (measurements performed using DNA Nanorulers, excited at 475 nm, imaged with Plan-Apochromat 63x / 1.4 Oil DIC and Axiocam 820 mono; resolution is sample and SNR dependent, optical resolution is diffraction limited).
Optical sectioning strength	Optical sectioning strength (z) down to 460 nm at 50% contrast (measured using Argolight SIM slide V2 excited at 475 nm, with Plan-Apochromat 63x/1.4 Oil DIC and Axiocam 820 mono; resolution is sample and SNR dependent)

Dimensions (width × depth × height)

Apotome 3 slider for Axio Imager	Approx. 278 mm × 90 mm × 76 mm
Apotome 3 slider for Axio Observer	Approx. 295 mm × 90 mm × 78 mm
Apotome 3 slider for Axio Zoom.V16	Approx. 278 mm × 90 mm × 76 mm

Weight

Apotome 3 slider	Approx. 1.1 kg
------------------	----------------

Functional Data

Area of use	Closed rooms
Radio interference suppression	As per EN 55011 Class A
Noise immunity	As per DIN EN 61326-1

Operating Data

Interference Suppression	In accordance with EN 55011 class A
Interference Resistance	In accordance with DIN EN 61326-1
Supply Voltage	24 V DC
Power Consumption Apotome 3	Max. 5 W

Technical Specifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › **Technology and Details**
- › Service

Ambient Conditions for Operation

Permissible ambient temperature	+5 to +30 °C
Permissible relative humidity	Max. 80 % at +30 °C
Air pressure	800 hPa to 1060 hPa
Operating altitude	Max. 2000 m
Pollution degree	2
Warm-up period	30 min

Grid Frequencies (transmission grid high / medium / low)

Axio Imager slider	5 / 9 / 17.5 lp / mm
Axio Observer slider	10 / 17.5 / 35 lp / mm
Axio Zoom.V16 slider	10 / 15 / 20 lp / mm

Installation Conditions

The grid projection method used for the Apotome 3 is sensitive to vibration, which can have various causes (including strong draughts). Vibrations are visible as streak artefacts in the resulting image. The microscope must therefore be set up so that it is exposed to as little vibration as possible on a vibration-damped table or on a suitable microscope base.



ZEISS Service – Your Partner at All Times

Your microscope system from ZEISS is one of your most important tools. For over 175 years, the ZEISS brand and our experience have stood for reliable equipment with a long life in the field of microscopy. You can count on superior service and support - before and after installation. Our skilled ZEISS service team makes sure that your microscope is always ready for use.

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › Technology and Details
- › **Service**

Procurement

- Lab Planning & Construction Site Management
- Site Inspection & Environmental Analysis
- GMP-Qualification IQ/OQ
- Installation & Handover
- IT Integration Support
- Startup Training

Operation

- Predictive Service Remote Monitoring
- Inspection & Preventive Maintenance
 - Software Maintenance Agreements
 - Operation & Application Training
 - Expert Phone & Remote Support
 - Protect Service Agreements
 - Metrological Calibration
 - Instrument Relocation
 - Consumables
 - Repairs

New Investment

- Decommissioning
- Trade In

Retrofit

- Customized Engineering
 - Upgrades & Modernization
- Customized Workflows via ZEISS arivis Cloud



Please note: Availability of services depends on product line and location



Carl Zeiss Microscopy GmbH
07745 Jena, Germany
microscopy@zeiss.com
www.zeiss.com/apotome

Follow us on social media:

