

# ZEISS VISTHESIA

## ZEISS VISTHESIA intra

### Technical Specifications



Seeing beyond

	VISTHESIA® 1.0 % <sup>1</sup>	VISTHESIA 1.5 % <sup>1</sup>
	VISTHESIA 1.0 % intra <sup>2</sup>	VISTHESIA 1.5 % intra <sup>2</sup>
Behavior	Viscous-cohesive	Viscous-cohesive
Origin	Bacterial fermentation	Bacterial fermentation
Substance concentration I (%)	1.0 Sodium Hyaluronate 1.0 Lidocaine Hydrochloride	1.5 Sodium Hyaluronate 1.0 Lidocaine Hydrochloride
Substance concentration II (mg/ml)	10 Sodium Hyaluronate 10 Lidocaine Hydrochloride	15 Sodium Hyaluronate 10 Lidocaine Hydrochloride
Volume (ml)	0.8	0.8
CDI	39	44
Osmolality (mOsmol/kg)	280–330	280–330
Molecular weight (megadalton <sup>4</sup> )	2.9	2.9
Pseudoplasticity	75	133
Zero-shear viscosity (Pa s <sup>5</sup> )	72	295
Storage temperature (°C)	2–8	2–8
Content VISTHESIA:	1 × intracameral syringe, 1 × 27 G cannula, 2 × topical ampules (see below)	1 × intracameral syringe, 1 × 27 G cannula, 2 × topical ampules (see below)
Content VISTHESIA intra:	1 × intracameral syringe, 1 × 27 G cannula	1 × intracameral syringe, 1 × 27 G cannula

	VISTHESIA topical ampule <sup>3</sup>
Substance concentration I (%)	2 × 0.3 Sodium Hyaluronate, 2.0 Lidocaine Hydrochloride
Substance concentration II (mg/ml)	2 × 3 Sodium Hyaluronate, 20 Lidocaine Hydrochloride
Volume (ml)	2 × 0.3
Osmolality (mOsmol/kg)	270–320
Molecular weight (megadalton <sup>4</sup> )	2.9
Zero-shear viscosity (Pa s <sup>5</sup> )	1
Storage temperature (°C)	2–8
Content:	2 × topical ampules

All OVDs are latex and phthalate free, have a shelf life of 3 years and a pH range between 6.8 and 7.6.


<sup>1</sup> VISTHESIA 1.0% and VISTHESIA 1.5% are not for sale in the UK or Portugal

<sup>2</sup> VISTHESIA intra version does not contain 2 topical ampules

<sup>3</sup> Not available individually; only available with VISTHESIA 1.0% or 1.5%

<sup>4</sup> Megadalton = 1 million Dalton

<sup>5</sup> Pa s = 1000 mPa s



**CE** 0344

VISTHESIA 1.0%  
VISTHESIA 1.5%  
VISTHESIA 1.0% intra  
VISTHESIA 1.5% intra



**Hyaltech Ltd.**

Starlaw Business Park  
Livingston EH54 8SF  
United Kingdom  
[www.zeiss.com/ovd](http://www.zeiss.com/ovd)  
[www.zeiss.com/contact](http://www.zeiss.com/contact)



**Carl Zeiss Meditec AG**

Goeschwitzer Str. 51–52  
07745 Jena  
Germany  
[www.zeiss.com/ovd](http://www.zeiss.com/ovd)  
[www.zeiss.com/contact](http://www.zeiss.com/contact)

**en-0US\_32\_022\_01391** Printed in Germany. CZ-IX/2021 International edition: Only for sale in selected countries.  
The contents of the brochure may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information. Subject to changes in design and scope of delivery and due to ongoing technical development. VISTHESIA 1.0%. VISTHESIA 1.5%. VISTHESIA 1.0% intra and VISTHESIA 1.5% intra are registered trademarks of Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and/or other countries.  
© Carl Zeiss Meditec AG, 2021. All rights reserved.