

# ZEISS AT TORBI 719M

# ZEISS AT TORBI 719MP

## Technical Specifications



Seeing beyond



### AT TORBI® 719M

Optic Design	Monofocal, bitoric, aspheric (aberration neutral)	
Material	Hydrophilic acrylic (25 %) with hydrophobic surface properties	
Optic Diameter	6.0 mm	
Total Diameter	11.0 mm	
Haptic Angulation	0°	
Lens Design	Single-piece, bitoric, MICS	
Incision Size	1.8 mm from -8.0 to +28.0 D Spherical Equivalent (SE) with +1.0 to +4.0 D cylinder / 2.2 mm for all remaining diopter-cylinder-combinations up to +32.0 D	
Company Labeled A-Constant <sup>1</sup>	118.4	
Diopter Range	Spherical Equivalent (SE) Cylinder	-4.0 to +32.0 D, 0.5 D increments <sup>2</sup> +1.0 to +12.0 D, 0.5 D increments <sup>2</sup>
ACD <sup>1</sup>	5.14	
Abbe Number	58-60	
Refractive Index	1.46	
Implantation in	Capsular Bag	
Injector/Cartridge Set	For IOLs from -8.0 to +28.0 D SE with +1.0 to +4.0 D cylinder*: AT.Shooter™ A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge Set or VISCOJECT-BIO 1.8 Injector Set	
	* For all remaining diopter-cylinder-combinations use VISCOJECT-BIO 2.2 Injector Set or AT.Shooter A1-2000 injector with VISCOJECT-BIO 2.2 Cartridge Set.	



### AT TORBI 719MP

Optic Design	Monofocal, bitoric, aspheric (aberration neutral)	
Material	Hydrophilic acrylic (25 %) with hydrophobic surface properties	
Optic Diameter	6.0 mm	
Total Diameter	11.0 mm	
Haptic Angulation	0°	
Lens Design	Single-piece, bitoric, MICS	
Incision Size	1.8 mm	
Company Labeled A-Constant <sup>1</sup>	118.4	
Diopter Range	Spherical Equivalent (SE) Cylinder	-8.0 to +28.0 D, 0.5 D increments <sup>2</sup> +1.0 to +4.0 D, 0.5 D increments <sup>2</sup>
ACD <sup>1</sup>	5.14	
Abbe Number	58-60	
Refractive Index	1.46	
Implantation in	Capsular Bag	
Injector/Cartridge Set <sup>3</sup>	BLUEMIXS® 180	

<sup>1</sup> For optimized A Constants and ACD Constants refer to IOLCon: [www.iolcon.org](http://www.iolcon.org).

<sup>2</sup> Further preselected SE/cylinder combinations are available above and below the stated SE range.

<sup>3</sup> Please refer to our web pages for the most up-to-date references.



0297

AT TORBI 719M  
AT TORBI 719MP  
BLUEMIXS 180 Injector  
AT.Shooter A1-2000



**Carl Zeiss Meditec AG**  
Goeschwitzer Strasse 51–52  
07745 Jena  
Germany  
[www.zeiss.com/iol](http://www.zeiss.com/iol)  
[www.zeiss.com/toric-iols](http://www.zeiss.com/toric-iols)



0482

VISCOJECT



**Medicel AG**  
Dornierstrasse 11  
9423 Altenrhein  
Switzerland  
[www.zeiss.com/injectors](http://www.zeiss.com/injectors)  
[www.zeiss.com/med/contacts](http://www.zeiss.com/med/contacts)

**en-OUS\_32\_022\_015111** Printed in Germany. CZ-XII/2022 International edition: Only for sale in selected countries.  
The contents of the brochure may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information.  
Subject to changes in design and scope of delivery and due to ongoing technical development. AT TORBI, AT.Shooter and BLUEMIXS are either trademarks or registered trademarks of  
Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and /or other countries. VISCOJECT is a trademark of Medicel (CE 0482).  
© Carl Zeiss Meditec AG, 2022. All rights reserved.