



Visalis 100 from ZEISS
Experience the difference





Test it!

Before making an important decision, it is always a good idea to test your options first. After all, nothing beats firsthand experience. That's important whether you're looking for a new car ...

Visalis 100 from ZEISS – Release 2.0

Experience phaco efficiency

... or for a new phaco machine. After experiencing the release 2.0 of the Visalis® 100 from ZEISS, surgeons provided exciting feedback on phaco efficiency:

95 %

of cases showed outstanding followability*

100 %

experienced superior chamber stability*

98 %

of hard cataract cases confirmed extremely good ultrasound performance*

The ZEISS Visalis 100 is unique in its class and provides everything you need and more than you would expect! ZEISS quality in phaco technology that you simply have to test to ...

... experience the difference.



* Data on file

Experience unique features for your patients ...



Capturing followability*

Its capturing followability attracts and continuously engages nucleus fragments on the phaco tip.

Balanced stroke – Adaptive Power Control

Microprocessor controls, measures and stabilizes the tip stroke intraoperatively, resulting in potentially less stress for your patient's eye.

Surge Security System

Intended to maintain chamber stability at occlusion break thereby minimizing risk of surge.

Minimize energy – new pulse modes

New pulse modes allow you adjust the on time according to your needs and thereby minimize the energy applied to the patient's eye.

Minimal induced astigmatism – CO-MICS capabilities

Unique in its class, the ZEISS Visalis 100 enables a 1.8 mm incision. As described in the literature, MICS reduces surgically induced astigmatism, increases patient comfort and accelerates wound healing.

Unique pneumatic guillotine cutter

The ZEISS Visalis 100 is the only system in its class with a pneumatic guillotine cutter that includes an integrated quiet compressor – a unique feature sharply focused on your patient's eye.

* Data on file

... and usability and cost-effectiveness for you!



Total control – aspiration rise time setting

Keep control of aspiration and vacuum with the new peristaltic pump and aspiration rise time setting that provides a venturi-like effect.

Low running costs – durability, reusability

ZEISS stands for high quality products that are engineered for long-term use. This, together with the wide variety of reusable accessories, ensures that the ZEISS Visalis 100 is tailored to meet the economic needs of your practice.



Excellent usability –

user interface, footswitch and prime options

Intuitive handling; new graphical user interface; new prime options; and a new, more ergonomic footswitch with simultaneous linear control of aspiration and ultrasound. A difference in usability you'll experience every day!

A winning team – combine with OPMI LUMERA 300

Combine the Visalis 100 with the OPMI LUMERA® 300 surgical microscope from ZEISS and get proven ZEISS optics with a clear, brilliant and high-contrast image together with a unique phaco system – a winning team for cataract surgery.

*ZEISS quality in cataract surgery –
ZEISS OPMI LUMERA 300 best-in-class visualization
meets the phaco efficiency of ZEISS Visalis 100*

Technical Data

Visalis 100 from ZEISS

Dimensions	Height: 205 mm, width: 400 mm, depth: 405 mm
Weight	12.5 kg
Irrigation	Gravity fed fluid delivery; eye pressure determined by height of irrigation source Solenoid valving element
Aspiration	Peristaltic pump 5 to 500 mmHg vacuum range Flow rate from 2 to 50 cc/min Surgery mode with fixed or linear control from 0 to preset of the flow and/or vacuum level via Double Linear Footswitch II 0.5 to 12 s aspiration rise time
Phacoemulsification	4 Piezoelectric crystals handpiece with frequency approx. 40 kHz U/S mode with Continuous, Pulse, Burst and with fixed or linear control from 0 to preset U/S Timer /Timer EPT from 00:00:0 minutes to 59:59:9 minutes
Vitreotomy	Handpieces type with pneumatically powered guillotine cutter Cutting rate from 60 to 700 cuts per minute Actuating medium with pressurized air from integrated air compressor
Diathermy	Adjustable bipolar power from 9 W @ 200 Ω , 5 to 100 % Linear control with Double Linear Footswitch II

Experiencing the difference of superior phaco efficiency.

ZEISS Visalis 100



// INNOVATION
MADE BY ZEISS

CE 0297



Carl Zeiss Meditec AG
Goeschwitzer Strasse 51–52
07745 Jena
Germany
www.zeiss.com/visalis100

EN_32_010_0030III Printed in Germany. CZ-VI/2017 International version.
Not all products, services or offers are approved or offered in every market and approved labeling and instructions may vary from one country to another. The contents of the brochure may differ from the current status of approval of the product in your country. Please contact our regional representative for more information. Subject to changes in design and scope of delivery and due to ongoing technical development. VISALIS and OPWII LUMERA are either trademarks or registered trademarks of Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and/ or other countries.
© Carl Zeiss Meditec AG, 2017. All rights reserved.