



ZEISS OPMI pico

Product specifications

Last updated: 2020-05



Seeing beyond

Product specifications

OPMI® pico from ZEISS is a compact, high-performance, easy-to-use surgical microscope:

- High-intensity fiber optic illumination
- Ergonomic arrangement of all controls
- Objective lens with fine focusing
- 5-step magnification changer
- Widefield optics
- Very mobile floor stand with four easy-to-move, lockable casters
- Suitable for SAD (Sexual Assault Documentation)



ZEISS OPMI pico is available in various configurations:

- Straight or inclined binocular tube for excellent ergonomics
- You can choose between three integrated HD video camera versions:
HD-ready 720p, Full HD 1080p or live & streaming (Full HD 1080p) with recording license.

Technical data

S100 / OPMI pico from ZEISS

| | | |
|--|---|-----------------|
| Magnification | Manual magnification changer with five settings | |
| Eyepieces | 12.5x widefield eyepieces | without reticle |
| | | with reticle |
| | 10x widefield eyepieces | without reticle |
| | | with reticle |
| Viewing tube | Straight binocular tube, f=170 mm <i>or</i> 45° inclined binocular tube <i>or</i> 180° swiveling binocular tube | |
| Objective lens | 300 mm, focusable <i>or</i> 250 mm, focusable <i>or</i> Varioskop® 100, focusable working distance 200–300 mm | |
| Illumination | Halogen or LED | |
| Red-free filter | Swing-in | |
| Integrated HD video camera (option) | HD-ready 720p <i>or</i> Full HD 1080p <i>or</i> Live & streaming (Full HD 1080p) with recording license | |
| S100 suspension mount | Floor stand <i>or</i> Wall mount <i>or</i> Ceiling mount | |



S100 / OPMI pico



Carl Zeiss Meditec AG
Goeschwitzer Strasse 51–52
07745 Jena
Germany
www.zeiss.com/gyn
www.zeiss.com/med/contacts