



**VISALIS 500 Family and  
Visalis 100 from ZEISS**  
Accessory catalog






Seeing beyond

# Table of contents

Page





I/A Tubing	3
Cataract Accessories	4
<b>Phacoemulsification</b>	<b>4</b>
Handpiece	4
Phaco Sets	5
Tips	7
Additional Accessories	7
Sleeves	8
<b>Irrigation/Aspiration</b>	<b>8</b>
Bimanual Technique	8
Coaxial Technique	9
<b>Diathermy</b>	<b>11</b>
<b>Anterior Vitrectomy</b>	<b>11</b>
Vitreoretinal Accessories	12
<b>Trocar and Infusion Systems</b>	<b>12</b>
<b>Posterior Vitrectomy</b>	<b>12</b>
<b>Endoillumination</b>	<b>13</b>
<b>Injection and Removal Systems</b>	<b>14</b>
<b>Endodiathermy</b>	<b>15</b>
Miscellaneous Accessories	16
System Options	17

## I/A Tubing VISALIS 500 Family\*

single-use		<p><b>I/A TUBING CASSETTE QUICKSET (40 pcs)</b></p> <p>Cassette I/A tubing "QuickSet" with drainage bag</p> <p>303060-0204-000</p>
single-use		<p><b>CATARACT ADVANCED IRRIGATION TUBING (10 pcs)</b></p> <p>Directs balanced salt solution from irrigation source to the surgical handpieces at predefined pressure during eye surgery on anterior segment (10 pcs)</p> <p>303060-0208-000</p>
single-use		<p><b>RETINA ADVANCED IRRIGATION TUBING (10 pcs)</b></p> <p>Directs balanced salt solution from irrigation source to the surgical handpieces at predefined pressure during eye surgery on posterior segment (10 pcs)</p> <p>303060-0210-000</p>

\*VISALIS® 500 Family includes VISALIS S500 and VISALIS V500

## I/A Tubing Visalis 100

reusable		<b>I/A TUBING SET, REUSABLE</b>	303060-0200-000
single-use		<b>DRAINAGE BAG</b> (10 pcs) <sup>1</sup> Non Sterile, single-use	303060-0203-000
single-use		<b>I/A TUBING SET</b> (10 pcs) single-use	303060-0201-000
single-use		<b>IRRIGATION ADMINISTRATION SET</b> (10 pcs)	303060-0205-000

## Cataract Accessories VISALIS 500 Family and Visalis 100

### Phacoemulsification

#### Handpiece







reusable		<b>ULITE PHACO HANDPIECE</b> Phaco Handpiece 4 piezoelectric, 40 KHz Tip stroke from 0 to 100 µm	303060-0001-000
----------	---	---	-----------------



# Cataract Accessories VISALIS 500 Family and Visalis 100

## Phacoemulsification

Phaco Sets (not available in all countries)







single-use		<b>PHACO SET PLUS 20G, Straight, 15° (10 pcs)</b>	20G Phaco Tip Straight, 15° 20G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 2.6-2.8	303060-1020-000
single-use		<b>PHACO SET PLUS 20G, Straight, 30° (10 pcs)</b>	20G Phaco Tip Straight, 30° 20G Silicone Sleeve (2 pcs) Test Chamber and Tip Wrench	RIS* 2.6-2.8	303060-1021-000
single-use		<b>PHACO SET PLUS 20G, Flared, 30° (10 pcs)</b>	20G Phaco Tip Flared 30° 20G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 2.6-2.8	303060-1022-000
single-use		<b>PHACO SET PLUS 21G, Biconical, 15° (10 pcs)</b>	21G Phaco Tip Biconical 15° 21G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 2.2-2.5	303060-1023-000
single-use		<b>PHACO SET PLUS 21G, Biconical, 30° (10 pcs)</b>	21G Phaco Tip Biconical 30° 21G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 2.2-2.5	303060-1024-000
single-use		<b>PHACO SET PLUS 21G, Flared, 30° (10 pcs)</b>	21G Phaco Tip Flared 30° 21G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 2.2-2.5	303060-1025-000

\* Recommended incision size in mm

# Cataract Accessories VISALIS 500 Family and Visalis 100

## Phacoemulsification

### Phaco Sets (not available in all countries)









single-use		<b>PHACO SET PLUS 22G, Biconical, 15° (10 pcs)</b>	22G Phaco Tip Biconical 15° 22G Silicone Sleeve (2 pcs) Test Chamber and Tip Wrench	RIS* 1.7-2.1	303060-1026-000
single-use		<b>PHACO SET PLUS 22G, Biconical, 30° (10 pcs)</b>	22G Phaco Tip 30° Biconical 22G Silicone Sleeve (2 pcs) Test Chamber and Tip Wrench	RIS* 1.7-2.1	303060-1027-000
single-use		<b>PHACO SET PLUS 22G, Flared, 30° (10 pcs)</b>	22G Phaco Tip Flared 30° 22G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 1.7-2.1	303060-1028-000
single-use		<b>SLEEVE, TEST CHAMBER SET PLUS 20G (10 pcs)</b>	20G Silicone Sleeve (2 pcs) and Test Chamber	RIS* 2.6-2.8	303060-1030-000
single-use		<b>SLEEVE, TEST CHAMBER SET PLUS 21G (10 pcs)</b>	21G Silicone Sleeve (2 pcs) and Test Chamber	RIS* 2.2-2.5	303060-1031-000
single-use		<b>SLEEVE, TEST CHAMBER SET PLUS 22G (10 pcs)</b>	22G Silicone Sleeve (2 pcs) and Test Chamber	RIS* 1.7-2.1	303060-1032-000

\* Recommended incision size in mm

# Cataract Accessories VISALIS 500 Family and Visalis 100

## Phacoemulsification



### Tips

reusable		19G PHACO TIP, 15°	Straight 15°	303060-0050-000	RIS* 2.9-3.2	Sleeve 303060-0100-000
		19G PHACO TIP, 30°	Straight 30°	303060-0051-000		
		19G PHACO TIP, 45°	Straight 45°	303060-0052-000		
reusable		20G PHACO TIP, 15°	Straight 15°	303060-0054-000	RIS* 2.6-2.8	Sleeve 303060-0101-000
		20G PHACO TIP, 30°	Straight 30°	303060-0055-000		
		20G PHACO TIP, 45°	Straight 45°	303060-0056-000		
reusable		20G PHACO TIP, BICONICAL	Biconical 15°	303060-0053-000	RIS* 2.6-2.8	Sleeve 303060-0101-000
reusable		20G PHACO TIP, 20°, ANGLED	Angled 20°	303060-0062-000	RIS* 2.6-2.8	Sleeve 303060-0101-000
reusable		21G PHACO TIP, BICONICAL, 15°	Biconical 15°	303060-0059-000	RIS* 2.2-2.5	Sleeve 303060-0102-000
		21G PHACO TIP, BICONICAL, 30°	Biconical 30°	303060-0060-000		
reusable		22G PHACO TIP, BICONICAL, 15°	Biconical	303060-0063-000	RIS* 1.7-2.1	Sleeve 303060-0103-000
		22G PHACO TIP, BICONICAL, 30°	Biconical 30°	303060-0065-000		
reusable		22G PHACO TIP, 30°, SLEEVELESS BIMANUAL	Straight 30°	303060-0057-000	RIS* 1.0-1.2	
reusable		22G PROGR. PHACO TIP, SLEEVELESS BIMANUAL	Flared 30°	303060-0058-000	RIS* 1.0-1.2	

### Additional tip for VISALIS V500

reusable		PARS PLANA PHACO TIP	20G (Posterior Fragmentation)	303060-0064-000		
----------	---	----------------------	-------------------------------	-----------------	--	--

### Additional Accessories

reusable		TEST CHAMBER (2 pcs) <sup>1</sup>				303060-0202-000
reusable		PHACO TIP WRENCH <sup>1</sup>				303060-0121-000

\* Recommended incision size in mm

# Cataract Accessories VISALIS 500 Family and Visalis 100



## Phacoemulsification

### Sleeves

reusable		<b>19G SILICONE SLEEVE (GREY) 2.9–3.2 mm (3 pcs)</b>	RIS* 2.9-3.2	303060-0100-000
reusable		<b>20G SILICONE SLEEVE (BLUE) 2.6–2.8 mm (3 pcs)</b>	RIS* 2.6-2.8	303060-0101-000
reusable		<b>21G SILICONE SLEEVE (YELLOW) 2.2–2.5 mm (3 pcs)</b>	RIS* 2.2-2.5	303060-0102-000
reusable		<b>22G SILICONE SLEEVE (GREEN) 1.7–2.1 mm (3 pcs)</b>	RIS* 1.7-2.1	303060-0103-000

## Irrigation/Aspiration

### Bimanual Technique

single-use			<b>IRR &amp; ASP HANDPIECES BIMANUAL (10 pcs)</b> Irrigation Side Ports & Aspiration Bimanual Handpieces Aspiration: 21G, 0.35 mm Irrigation: 21G, 2 Side Ports 0.45 mm	303060-0165-000
single-use			<b>IRR &amp; ASP HANDPIECES BIMANUAL FRONT (10 pcs)</b> Irrigation Front Port & Aspiration Bimanual Handpieces Aspiration: 21G, 0.35 mm Irrigation: 21G, 1 Front Port 0.5 mm	303060-0166-000
reusable			<b>IRRIGATION HANDPIECE, BIMANUAL</b> 21G, 0.6 mm	303060-0155-000
reusable			<b>ASPIRATION HANDPIECE, BIMANUAL</b> 21G, 0.3 mm	303060-0156-000

\* Recommended incision size in mm

# Cataract Accessories VISALIS 500 Family and Visalis 100

## Irrigation/Aspiration

### Coaxial Technique














single-use		<b>I/A HANDPIECE, 20G STRAIGHT</b> 20G I/A Handpiece (10 pcs), Coaxial, Straight	RIS* 2.6-2.8	303060-0141-000
single-use		<b>I/A HANDPIECE, 20G ANGLED 45°</b> 20G I/A Handpiece (10 pcs), Coaxial, Angled 45°	RIS* 2.6-2.8	303060-0142-000
single-use		<b>I/A HANDPIECE, 21G ANGLED 45°</b> 21G I/A Handpiece (10 pcs), Coaxial, Angled 45°	RIS* 2.2-2.5	303060-0143-000
single-use		<b>I/A HANDPIECE, 21G CURVED</b> 21G I/A Handpiece (10 pcs), Coaxial, Curved 40°	RIS* 2.2-2.5	303060-0146-000
single-use		<b>I/A HANDPIECE, 22G ANGLED 45°</b> 22G I/A Handpiece (10 pcs), Coaxial, Angled 45°	RIS* 1.7-2.1	303060-0144-000
single-use		<b>I/A HANDPIECE, 22G CURVED</b> 22G I/A Handpiece (10 pcs), Coaxial, Curved 40°	RIS* 1.7-2.1	303060-0145-000
reusable		<b>I/A HANDPIECE, COAXIAL, REUSABLE</b> I/A Handpiece, Coaxial		303060-0140-000

\* Recommended incision size in mm

# Cataract Accessories VISALIS 500 Family and Visalis 100

## Irrigation/Aspiration




### Coaxial Technique

reusable		I/A CANNULA, METAL SLEEVE, STRAIGHT, Ø 1.6 mm		303060-0150-000
		I/A CANNULA, METAL SLEEVE, ANGLED, Ø 1.6 mm		303060-0151-000
		I/A CANNULA, METAL SLEEVE, CURVED, Ø 1.6 mm		303060-0152-000
		I/A CANNULA, METAL SLEEVE, 12 O'CLOCK, Ø 1.6 mm		303060-0153-000
		I/A CANNULA, METAL SLEEVE, 90°, Ø 1.6 mm		303060-0154-000
reusable		19G I/A CANNULA, SILICONE SLEEVE, STRAIGHT	RIS* 2.9-3.2	303060-0157-000
		19G I/A CANNULA, SILICONE SLEEVE, ANGLED		303060-0158-000
reusable		20G I/A CANNULA, SILICONE SLEEVE, STRAIGHT	RIS* 2.6-2.8	303060-0159-000
		20G I/A CANNULA, SILICONE SLEEVE, ANGLED		303060-0160-000
reusable		21G I/A CANNULA, SILICONE SLEEVE, STRAIGHT	RIS* 2.2-2.5	303060-0161-000
		21G I/A CANNULA, SILICONE SLEEVE, ANGLED		303060-0162-000
reusable		22G I/A CANNULA, SILICONE SLEEVE, STRAIGHT	RIS* 1.7-2.1	303060-0163-000
		22G I/A CANNULA, SILICONE SLEEVE, ANGLED		303060-0164-000

\* Recommended incision size in mm

## Cataract Accessories VISALIS 500 Family and Visalis 100

### Diathermy

reusable		<b>DIATHERMY PENCIL ERASER</b>	303060-0252-000
reusable		<b>DIATHERMY FORCEPS</b>	303060-0251-000
reusable		<b>DIATHERMY BIPOLAR CABLE<sup>1</sup></b>	303060-0250-000

### Anterior Vitrectomy

single-use		<b>20G A-VIT PROBE WITH SLEEVE (5 pcs)</b> 20G Anterior Vitrectomy Probe Silicone Sleeve	303060-0301-000
------------	---	--	-----------------









# Vitreoretinal Accessories VISALIS V500

## Trocar and Infusion Systems

single-use		RETILOCK® 23G NON VALVED SELF-RETAINING TROCAR (5 pcs)	303060-0414-000*
single-use		RETILOCK® 23G VALVED SELF-RETAINING TROCAR (5 pcs)	303060-0415-000*
single-use		RETILOCK® 25G NON VALVED SELF-RETAINING TROCAR (5 pcs)	303060-0416-000*
single-use		RETILOCK® 25G VALVED SELF-RETAINING TROCAR (5 pcs)	303060-0417-000*
single-use		20G SCLERAL INFUSION CANNULA (5 pcs)	303060-0503-000
single-use		23G TROCAR KIT ONE-STEP (10 pcs)	303060-0405-000**
single-use		23G TROCAR KIT ONE-STEP OIL INFUSION (10 pcs)	303060-0406-000**
single-use		25G TROCAR KIT ONE-STEP (10 pcs)	303060-0413-000**
single-use		25G TROCAR KIT ONE-STEP OIL INFUSION (10 pcs)	303060-0412-000**
single-use		25G TROCAR KIT (5 pcs)	303060-0409-000**

\* For sale in Europe only. \*\* No longer available in Europe.











## Posterior Vitrectomy

single-use		 20G POSTERIOR VITRECTOMY PROBE (5 pcs) High Speed	303060-0403-000
single-use		 23G POSTERIOR VITRECTOMY PROBE (5 pcs) High Speed	303060-0407-000
single-use		 25G POSTERIOR VITRECTOMY PROBE (5 pcs) High Speed	303060-0408-000








# Vitreoretinal Accessories VISALIS V500

## Endoillumination

single-use		 <b>20G ENDO-ILLUMINATION PROBE (5 pcs)</b> 303060-0600-000
		 <b>20G WIDE-ANGLE ENDO-ILLUMINATION PROBE (5 pcs)</b> 303060-0601-000
single-use		 <b>20G ENDO-ILLUMINATION PROBE MULTIPOINT (5 pcs)</b> 303060-0602-000
		 <b>20G WIDE-ANGLE ENDO-ILLUMIN. MULTIPOINT (5 pcs)</b> 303060-0603-000
single-use		 <b>23G ENDO-ILLUMINATION PROBE MULTIPOINT (5 pcs)</b> 303060-0604-000
		 <b>23G WIDE-ANGLE ENDO-ILLUMIN. MULTIPOINT (5 pcs)</b> 303060-0605-000
single-use		<b>25G ENDO-ILLUMINATION PROBE MULTIPOINT (5 pcs)</b> 303060-0606-000





# Vitreoretinal Accessories VISALIS V500

## Injection and Removal Systems




single-use		TWO WAY STOPCOCK (5 pcs)	303060-0505-000
single-use		ADAPTER AIR/CONTROLLED IRRIGATION (5 pcs)	303060-0506-000
single-use		AIR INJECTION TUBE WITH FILTER (5 pcs)	303060-0508-000
single-use		SILICONE OIL INJECTION TUBE (5 pcs)	303060-0509-000
single-use		VISCIOUS FLUID REMOVAL KIT (10 pcs)	303060-0507-000

## Vitreoretinal Accessories VISALIS V500


### Endodiathermy

single-use		20G ENDODIATHERMY PROBE (5 pcs)	303060-0262-000
single-use		23G ENDODIATHERMY PROBE (5 pcs)	303060-0263-000
single-use		25G ENDODIATHERMY PROBE (5 pcs)	303060-0264-000
reusable		DIATHERMY BIPOLAR CABLE <sup>1</sup>	303060-0250-000


## Miscellaneous Accessories VISALIS 500 Family

single-use		<p><b>SCREEN DRAPE VISALIS 500</b> (8 pcs)</p>	303060-0352-000
single-use		<p><b>REMOTE CONTROL COVER VISALIS 500</b> (10 pcs)</p>	303060-0353-000
single-use		<p><b>TRAY COVER VISALIS 500</b> (2 pcs)</p>	303060-0354-000

## Miscellaneous Accessories VISALIS 500 Family and Visalis 100


reusable		<p><b>STERILIZATION TRAY<sup>1</sup></b> 26,5 x 16,3 x 4,6 cm</p>	303060-0350-000
----------	---	---	-----------------

## Miscellaneous Accessories Visalis 100

single-use		<p><b>SCREEN COVER &amp; MAYO DRAPE</b> (2 pcs) Sterile</p>	303060-0351-000
------------	---	---	-----------------

## System Options VISALIS 500 Family

---




**Console Tray<sup>2</sup>**  
with manually adjustable height

303041-9004-000

---

## Additional System Option for VISALIS V500

---



**TWO-WAY ADAPTER XENON-LIGHTSOURCE<sup>2</sup>**

303060-0608-000

---



All products without  
superscript number



All products with the  
superscript number 1

**Registered with VISALIS 500 Family**  
All products with the superscript number 2



**Carl Zeiss Meditec AG**  
Goeschwitzer Strasse 51–52  
07745 Jena  
Germany  
[www.zeiss.com/med/contacts](http://www.zeiss.com/med/contacts)  
[www.zeiss.com/visalis500](http://www.zeiss.com/visalis500)

**EN\_32\_010\_0113VI** Printed in Germany

The contents of the brochure may differ from the current status of approval of the product in your country. Please contact our regional representative for more information. Subject to change in design and scope of delivery and as a result of ongoing technical development. VISALIS is a registered trademark of Carl Zeiss Meditec AG. RETILOCK is a registered trademark of FCI S.A.S (CE 0459). © Carl Zeiss Meditec AG, 2019. All copyrights reserved.