

ZEISS CT SPHERIS 209M

Technical Specifications



Seeing beyond



CT SPHERIS® 209M

Optic Design	Monofocal, spheric
Material	Hydrophilic acrylic (25%) with hydrophobic surface properties
Optic Diameter	6.0 mm
Total Diameter	11.0 mm
Haptic Angulation	0°
Lens Design	Single-piece, MICS
Incision Size	1.8 mm
Company Labeled A-Constant¹	117.9
Diopter Range	From 0.0 to +32.0 D 0.0 to +10.0 D, 1.0 D increments +10.0 to +30.0 D, 0.5 D increments +30.0 to +32.0 D, 1.0 D increments
ACD	4.82
Implantation in	Capsular Bag
Injector/Cartridge Set²	AT.Shooter A1-2000/VISCOJECT™-BIO 1.8 Cartridge Set <i>or</i> VISCOJECT-BIO 1.8 Injector Set <i>or</i> A6 Injector set

¹ For optimized A Constants and ACD Constants refer to IOLCon: www.iolcon.org

² Please refer to our web page for the most up-to-date references: www.zeiss.com/injectors



0297

AT.Shooter A1-2000
CT SPHERIS 209M



Carl Zeiss Meditec AG

Goeschwitzer Strasse 51–52
07745 Jena
Germany
www.zeiss.com/spheris
www.zeiss.com/med/contacts



0482

VISCOJECT-BIO 1.8 Cartridge Set
VISCOJECT-BIO 1.8 Injector Set



Medicel AG

Dornierstrasse 11
9423 Altenrhein
Switzerland
www.zeiss.com/injectors
www.zeiss.com/med/contacts



0120

A6 Injector Set



MDJ SAS

ZI Racine Palladuc
63650 La Monnerie Le Montel
France

en-OUS_32_022_016911 CZ-VIII/2023 International edition: Only for sale in selected countries.

The contents of the brochure may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information. Subject to changes in design and scope of delivery and due to ongoing technical development. AT.Shooter A1-2000 and CT SPHERIS are either trademarks or registered trademarks of Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and/or other countries. VISCOJECT is a trademark of Medicel (0482). A6 injector is a trademark of MDJ (0120). © Carl Zeiss Meditec AG, 2023. All rights reserved.