



CT ASPHINA 404
Technical Specifications



CT ASPHINA 404



CT ASPHINA® 404

| | |
|---|---|
| Optic Design | Monofocal, aspheric (aberration neutral) |
| Material | Hydrophilic acrylic (25 %) with hydrophobic surface properties |
| Optic Diameter | 6.0 mm |
| Total Diameter | 11.0 mm |
| Haptic Angulation | 0° |
| Lens Design | Single-piece |
| Incision Size | 2.2 mm |
| Company Labeled A-Constant ¹ | 118.3 |
| Diopter Range | From -10.0 to +30.0 D (from +31.0 to +42.0 D upon request) -10.0 to +10.0 D: 1.0 D increments +10.0 to +30.0 D: 0.5 D increments +30.0 to +42.0 D: 1.0 D increments |
| ACD | 5.14 |
| Implantation in | Bag |
| Injector/Cartridge Set ² | For IOLs from -10.0 to +30.0 D*: AT.Shooter A1-2000 / VISCOJECT™-BIO 2.2 Cartridge Set <i>or</i> VISCOJECT-BIO 2.2 Injector Set * For IOLs from +31.0 to +42.0 D ACCUJECT™ 3.0 Injector Set |

¹ Please refer to our web pages for optimized A-Constants.

² Please refer to our web pages for the most up-to-date references.



CT ASPHINA 404



Carl Zeiss Meditec AG
Goeschwitzer Str. 51–52
07745 Jena
Germany
www.zeiss.com/iol
www.zeiss.com/med/contacts