

# ZEISS COMBIVISC

## Technical Specifications



Seeing beyond

	<b>Cohesive part</b>	<b>Dispersive part</b>
Behavior	Viscous-cohesive	Medium-viscosity dispersive
Origin	Bacterial fermentation	Bacterial fermentation
Substance	Sodium Hyaluronate	Sodium Hyaluronate
Concentration I (%)	1.5	3.0
Concentration II (mg/ml)	15	30
Volume (ml)	1.0	0.85
CDI	44	18
Osmolality (mOsmol/kg)	300 – 360	300 – 360
Molecular weight (megadalton <sup>1</sup> )	2.9	1.0
Pseudoplasticity	133	10
Zero-shear viscosity (Pa s <sup>2</sup> )	295	47
Cannula	25 G	25 G
Storage temperature (°C)	2 – 8	2 – 8
Singlepack content:	2 × intracameral syringes, 2 × 25 G cannulas	
Multipack content:	10 × intracameral syringes, 10 × 25 G cannulas	

COMBIVISC® is free of natural dry latex and phthalate, has a shelf life of 3 years and a pH range between 6.8 and 7.6.

<sup>1</sup> Megadalton = 1 million Dalton

<sup>2</sup> Pa s = 1000 mPa s





0344

COMBIVISC



**Hyaltech Ltd.**

Starlaw Business Park  
Livingston EH54 8SF  
United Kingdom  
[www.zeiss.com/ovd](http://www.zeiss.com/ovd)  
[www.zeiss.com/contacts](http://www.zeiss.com/contacts)



**Carl Zeiss Meditec AG**

Goeschwitzer Str. 51–52  
07745 Jena  
Germany  
[www.zeiss.com/ovd](http://www.zeiss.com/ovd)  
[www.zeiss.com/contacts](http://www.zeiss.com/contacts)

**en-OUS\_32\_022\_0137111** Printed in Germany. CZ-III/2023 International edition: Only for sale in selected countries.  
The contents of the brochure may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information. Subject to changes in design and scope of delivery and due to ongoing technical development.  
COMBIVISC is either trademark or registered trademark of Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and/or other countries.  
© Carl Zeiss Meditec AG, 2023. All rights reserved.