

ZEISS CT LUCIA 202

Technical Specifications



Seeing beyond



CT LUCIA® 202

Optic Design	Monofocal, spheric
Material	Hydrophobic acrylic
Optic Diameter	6.0 mm
Total Diameter	13.0 mm
Haptic	5°
Lens Design	Three-piece
Incision Size	3 mm
Company Labeled A-Constant¹	118.2
Diopter Range	From +4.0 to +34.0 D, 0.5 D increments
ACD²	5.1
Abbe number	50
Implantation in	Capsular bag
Injector/Cartridge Set²	ACCUJECT 3.0 3P for diopter range +4.0 to +27.0 D for higher powers use forceps

¹ For optimized A Constants and ACD Constants refer to IOLCon: www.iolcon.org.

² Please refer to our web page for the most up-to-date references: www.zeiss.com/injectors

CE 0297

CT LUCIA 202



Carl Zeiss Meditec AG

Goeschwitzer Strasse 51–52

07745 Jena

Germany

www.zeiss.com/lucia

www.zeiss.com/med/contacts

en-INT_32_022_01931 CZ-III/2026 International edition: Only for sale in selected countries.
The contents of the datasheet may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information.
Subject to changes in design and scope of delivery and due to ongoing technical development. CT LUCIA is either a trademark or registered trademark of Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and/or other countries.
© Carl Zeiss Meditec AG, 2026. All rights reserved.