

ZEISS QUATERA 700

Unleashing a new experience in phaco



New Features

Elevating the ZEISS QUATERA 700 Experience

ZEISS 21G Angled Phaco Tip¹

- 21G single use tip
- Angled 20°
- Bevel 30°



There's more with the new SW 1.2.0 update

- Wireless foot control panel
- POD override
- New service GUI
- Sound slider (separately adjustable sound volumes)
- CALLISTO eye 5.0 support



For more information visit zeiss.com/quatera700

en-INT_32_010_01781 CZ-X/2024 International edition: Only for sale in selected countries.
The contents of this document may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information. Subject to changes in design and scope of delivery and due to ongoing technical development. QUATERA 700 is either a trademark or registered trademark of Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and / or other countries. © Carl Zeiss Meditec AG, 2024. All rights reserved.

¹ Material number -303060-1019-000



Seeing beyond