

Single-use Lenses for ZEISS RESIGHT

Uncompromised
intraoperative
retina visualization



Seeing beyond



The Single-use Lenses for the RESIGHT® fundus viewing system from ZEISS provide retinal surgeons with a clear and uncompromised view for every surgical procedure in the posterior segment, while also offering all the advantages of disposable lenses. You will benefit from high-resolution visualization of the far periphery, with high light transmission enabling you to operate under reduced illumination intensity.



zeiss.com/su-lenses-resight

Smooth workflow and compatibility

All Single-use Lenses for ZEISS RESIGHT are compatible with any ZEISS ophthalmic microscope configured with ZEISS RESIGHT 500 or ZEISS RESIGHT 700, as well as with ZEISS integrated intraoperative OCT.

All lenses are delivered in separate sterile packaging, allowing your OR staff to prepare independently for surgery.



Clear and detailed view for safe macula surgery with the MACULA LENS



- High resolution with an extensive field of view and enhanced contrast
- Better view of surgical instruments during the procedure

Visualization up to the ora serrata with the ULTRA WIDE-ANGLE LENS



- Reduced need for tilting and scleral indentation
- Unobstructed access to the surgical instruments and reduced condensation

Versatility in peripheral surgery with the WIDE-ANGLE LENS



- Adaptability across posterior segment procedures without the need for switching lenses
- Outstanding peripheral visualization in high resolution

Technical Specifications			
	MACULA LENS	WIDE-ANGLE LENS	ULTRA WIDE-ANGLE LENS
Material	Optics: Polymethylmethacrylat (PMMA); Lens holder: Acrylnitril-Butadien-Styrol Copolymer (ABS)		
Sterilized with	Ethylene oxide		
Diameter	~ 23 mm	~ 19 mm	~ 21-26 mm



RESIGHT 500
RESIGHT 700



0297

Single-use lenses for
RESIGHT



Carl Zeiss Meditec AG
Goeschwitzer Strasse 51–52
07745 Jena
Deutschland
www.zeiss.com/med/contacts

en-INT_32_010_01741 CZ-VIII/2024 International edition: Only for sale in selected countries.
The contents of the flyer may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information. Subject to changes in design and scope of delivery and due to ongoing technical development. RESIGHT is either a trademark or registered trademark of Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and / or other countries.
© Carl Zeiss Meditec AG, 2024. All rights reserved.