


Operate more efficiently with full control over your instrument

ZEISS SmartSEM

SmartSEM is an operating system for electron microscopes that provides access to advanced microscope settings, designed to solve even the most challenging tasks. Through its comprehensive toolbar and configuration options, the SmartSEM software allows you to monitor and operate most of the active components of the microscope. Thus, it builds the bridge between control software and microscope hardware. The SmartSEM program suite includes several programs and utilities to fit your individual needs.

Furthermore, the software offers the possibility to memorize user-specific preferences, which eliminates the need to reconfigure the user interface each time the SmartSEM is started.



Highlights

Enhanced imaging conditions

Benefit from various imaging settings to receive excellent results and perform imaging tasks up to your needs

Additional software modules

Profit from a wide range of optional software modules which provide additional application-specific functionalities

More information on your images

Receive topographical and compositional information by overlaying images from different detectors

Customized usage

Personalize the user interface through the fully customizable toolbar that can be modified to fit your needs

Upgrade your ZEISS microscope

Upgrade your ZEISS microscope with the latest SmartSEM Software and benefit from advanced microscope settings and additional software modules.

As technical requirements may apply on some systems, please contact us to learn more about SmartSEM and how your processes will benefit from an upgrade:


microscopy@zeiss.com



Seeing beyond

ZEISS SmartSEM

Operate more efficiently with your scanning electron microscope



The ZEISS SmartSEM user interface provides a variety of different settings for advanced microscopy.

Availability*

- Crossbeam series
- GeminiSEM series
- MERLIN series
- Sigma series
- EVO series
- AURIGA series
- ULTRA series
- SUPRA series
- NVision series
- NEON series
- 15xx XB series
- 15xx series
- 14xx series

*Please note that SmartSEM v6.05 or higher requires at least a 24" monitor

Not for therapeutic use, treatment or medical diagnostic evidence. Not all products are available in every country. Contact your local ZEISS representative for more information.
EN_43_024_159 | CZ-10/2020 | Design, scope of delivery and technical progress subject to change without notice. | © Carl Zeiss Microscopy GmbH



microscopy@zeiss.com
www.zeiss.com/microscopy