



Upgrade your ZEISS Workstation to Windows 11 X-Ray Microscopes



Seeing beyond



Enhance Security and Performance

zeiss.com/microscopy

Future-Proof Your ZEISS X-ray Microscope by upgrading to Windows 11

Benefit from enhanced data protection and optimized performance that accelerates your research.

Compatible ZEISS Workstations Ready for Windows 11

Upgrade on-site if you're using one of the following systems:

- Primary Workstation: HP Z8 G5
- Primary Workstation: HP Z8 G4
- High-Performance Workstation: HP Z8 G5
- High-Performance Workstation: HP Z8 G4

Why Upgrade to a New ZEISS Workstation?

If your current workstation isn't Windows 11 compatible, upgrading to the available new workstation brings powerful advantages:

- **24/7 Reliability:** Enjoy uninterrupted operation with robust SSDs and power supplies, designed for demanding microscopy tasks.
- **Proven Stability:** ZEISS workstations are extensively tested, validated, and delivered with frozen configurations for long-term stability and easy service.
- **Faster Processing:** Boost performance and reduce processing time with the latest high-speed hardware.
- **Enhanced Graphics Support:** Stay ahead with seamless compatibility for the latest ZEN navx, Advanced Reconstruction Toolbox (ART) software, and cutting-edge graphics cards.

Windows 11 and TPM**: A New Standard in Security

With the gradual end of Microsoft support for Windows 10, we strongly recommend upgrading your system to the current Windows operating system with full security support.

Windows 11 has the ability to incorporate Trusted Platform Module (TPM), which adds an extra layer of security by securely storing cryptographic keys and ensuring that only authorized software can access the system. This is particularly crucial for systems that control sensitive equipment or are used in specialized environments.

ZEISS Workstation and Software Upgrade: Future-Proofing and Enhanced Features

ZEISS offers a complete software upgrade package, including Windows 11, a powerful workstation and the latest version of ZEN navx, Scout-and-Scan, or ART for uncompromised user experience. With the latest software and hardware package, your XRM Microscope is ensured to operate at the highest level, accelerating your research.

ZEISS utilizes Long Term Servicing Channel (LTSC) Microsoft products, which are not available as standard commercial Windows options for general users. Unlike other Windows servicing models that regularly push feature updates, LTSC products remain stable and are replaced every few years with a newer, well-tested version.

For more information on how you can benefit from a full Windows 11 and workstation upgrade, please contact us at:



microscopy@zeiss.com
www.zeiss.com/microscopy

* Please contact your local ZEISS contact to check if additional upgrades are required to maintain your current VersaXRM functionality.

** Please contact your local ZEISS contact to check if TPM is available for your workstation.

Carl Zeiss Microscopy GmbH
07745 Jena, Germany

Email: microscopy@zeiss.com
www.zeiss.com/microscopy