

Accelerate Your Qualification Process


ZEISS IQ/OQ Service

The qualification of quality assurance systems is essential for complying with the requirements of regulated industries. Speed up the integration time of your new quality assurance system with our full-service package featuring the relevant documentation of the Installation Qualification (IQ) and Operational Qualification (OQ) for our hardware and software.

Your benefits: ✓

- Time saving – start production earlier
- Compliance – meet national and international requirements in the pharmaceutical or biotechnology industry
- Expertise – check completeness of relevant documents
- Expert consulting – benefit from a single data source and a designated contact

IQ and OQ are important phases within the validation process of new systems:



OQ: Manufacturers verify that equipment operates according to functional specifications. They need to carry out tests for each component to confirm that every part of the microscope operates as intended at pre-set thresholds.

IQ: Companies operating in regulated markets must verify that the microscope meets design specifications and has been correctly installed and configured according to requirements. Manufacturers must confirm that all required documents, certificates, and manuals are correctly executed and available.

Our IQ/OQ Service includes the requirements for validation and provides all necessary documents for internal validation and process approval:

- Personal consulting and support in the IQ/OQ process
- Collection of relevant IQ/OQ documents
- IQ** > Verification of the scope of delivery and installation, complete documentation, IQ report
- OQ** > System and functional tests, system calibration, OQ report

Looking to accelerate your validation process?

Contact our experts:

validate.your.microscopy@zeiss.com
www.zeiss.com/biopharma