


ZEISS Remote Health Check



Seeing beyond



Benefit from a scheduled check of your coordinate measuring machine

With the ZEISS Remote Health Check, you can proactively record the current condition of your coordinate measuring machine together with a ZEISS expert. The predictive approach of the ZEISS Remote Health Check provides you with specific recommendations for action regarding the technical condition of your coordinate measuring machine displayed to you.

From *service for the customer*, to *service with the customer* - working together to improve technical machine condition with our ZEISS Remote Health Check product.

Your benefits:

- Proactive analysis and diagnosis with a predictive recommendation from ZEISS experts.
- Securing measurement capability and increasing machine availability.
- Planned repairs instead of unforeseen issues
- Joint execution with our ZEISS experts
- Secure data connection - protected by ZEISS and controlled by you

Technical requirements:

- Temporary connection to your evaluation computer (preferably with TeamViewer).
- Operating system: Windows 10 with admin rights
- ZEISS Smart Services Health Check Kit available on site

probes.zeiss.com